Orthocryl® LC

CE 0483

- Gebrauchsanweisung | **Lichthärtender Kunststoff**
Instructions for use | **Light-curing acrylic**
Mode d'emploi | **Résine photopolymérisable**
Instrucciones de uso | **Acrílico fotopolimerizable**
Modalità d'uso | **Resina fotopolimerizzabile**

D
DENTAURUM
1886

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Dear customer,

Thank you for choosing a quality product from Dentaureum.

In order to achieve the best results with this product, it is important to carefully study and follow these Instructions for use.

Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaureum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaureum GmbH & Co. KG
Turnstraße 31 | 75228 Ispringen | Germany

2. General product description

Orthocryl® LC is a light-curing acrylic for removable and fixed orthodontic appliances. It can also be used to fabricate bite plates, anti-snoring devices and drilling templates for implants. Orthocryl® LC does not contain methyl methacrylate or dibenzoyl peroxide.

3. Form of presentation

Orthocryl® LC is delivered in opaque cartridges that fit the injector (REF 160-420-00) (Fig. 1). The Orthocryl® LC bite paste is available in opaque jars.

4. Preparation

It is mandatory to wear face protection/safety glasses (REF 162-920-00) and special protective gloves, e.g. Dermatril® disposable gloves made of special nitrile (REF 162-907-00, 162-908-00 and 162-909-00) to avoid direct skin contact when processing Orthocryl® LC. Please refer also to the corresponding safety data sheets.



- Remove the small blue protective cap on the Orthocryl® LC cartridge by turning it anti-clockwise (Fig. 2).
- Remove the large blue protective cap on the Orthocryl® LC cartridge by pushing up the protruding flap (Fig. 3).
- If necessary, shorten an interchangeable cartridge tip (REF 160-421-00) to the length required and rotate onto the top of the cartridge as far as it will go (Fig. 4).
- The button on the upper side of the injector releases the forward and back feed of the injector plunger. Hold the button down and bring the plunger into initial position.
- Put the cartridge onto the injector and lock with a 90° turn (Fig. 5).
- Push the plunger forward whilst holding the button down until resistance is felt.
- Apply the emerging Orthocryl® LC onto wires, screws and models with care.

4.1 Separate plaster model / plastic model

Both plaster and plastic models **must** be given a layer of separating medium before Orthocryl® LC is applied.

When working with plaster models, we recommend using a special separating medium (REF 162-800-00) in order to attain a smooth, silky basal surface on the appliance and to avoid discolorations and plaster inclusions. It is not necessary to soak the model beforehand.

Plastic models **must** be given a layer of separating medium before Orthocryl® LC is applied in order to prevent the plastics from sticking. We recommend the 3D separating medium (REF 162-801-00). This is, however, not available in the USA.

4.1.1 Application of the separating medium (REF 162-800-00) on plaster models

- Apply the separating medium from the dosing bottle (REF 162-100-00) directly on the model and spread evenly using a brush.
- Stand the model on its back (tuber plane) so that excess fluid can drain away.
- After approx. 30 seconds, carefully blow away separating medium residues or droplets with compressed air and, if necessary, remove residues from screws and wire retentions with a dry brush.
- Application of the acrylic can start immediately afterwards!

4.1.2 Application of the 3D separating medium (REF 162-801-00) on plastic models

To fabricate dental models, we recommend using high quality plastic models that were made using DLP (Digital Light Processing). e. g. prime4me® base model (REF 430-110-10). Besides being highly accurate and therefore ensuring very precise fitting, these models allow for the usual process in the laboratory.

Due to the large variety of materials for model printing, it is essential to test the suitability of the 3D separating medium before its use.

- Apply the 3D separating medium evenly with a brush over the entire model before the expansion screw (or screws) and wire elements are finally fixed in place and allow the separating medium to dry thoroughly. To this end, apply the separating medium from the dosing bottle (REF 162-100-00) directly on the model and spread evenly using a brush.
- Stand the model on its back so that excess fluid can drain away.
- Allow the layer of separating medium to dry thoroughly.
- The separating medium is completely dry when it no longer feels sticky.
- Repeat the process in order to achieve an even separating layer everywhere.
- In order to accelerate the drying of the separating medium or to shorten the drying phases, careful blowing with compressed air or a household hair dryer can be helpful. For faster drying of the separating medium, the model can also be placed in a light curing unit for a few minutes*.

* The time required in the light curing unit may vary depending on the size of the workpiece and the amount of separating medium. The drying time also depends on the type of device – with the Solidilite V unit from Shofu, up to 3 minutes should be enough.

- If required, block out undercuts in the tooth/jaw area and fix the expansion screw(s) and wire elements as usual on the dried plastic model.
- Application of the acrylic can start immediately afterwards!

5. Processing Orthocryl® LC

Example: expansion plate

Orthocryl® LC is generally applied to the model that is dry and isolated but not hydrated. Mark the position of the expansion screws on the model and drill a hole or a groove with a conical bur (e.g. REF 123-598-00) for the screw holder. Fill the hole on the model with heated wax (e.g. REF 119-150-00) and insert the expansion screw with the lower, basal part of the holder. Remove the expansion screw after the wax has hardened.

Fix bent wire elements to the model with sticky wax (REF 122-300-00).

Isolate the model as described under point 4.1.

Application of acrylic can start immediately afterwards!

Apply Orthocryl® LC to the removed expansion screw so that the screw is completely surrounded by the material (Fig. 6). Please observe that the expansion screw is completely embedded in acrylic without any air bubbles. Afterwards insert the expansion screw back into the drilled hole with the wax coating. The wax in the hole ensures that the expansion screw settles in the same position as determined before.

With wire retentions, apply Orthocryl® LC to the space between wire and plaster model first and then embed the complete wire (Fig. 7).

Next, apply Orthocryl® LC until the desired form and thickness of the appliance (maxilla approx. 3 mm / mandible approx. 4 – 5 mm) is achieved (Fig. 8). Should this process require more than a single step, dip the tip of the cartridge slightly into the material already applied before each new step. This prevents air inclusions between the layers.

Fig. 6


Fig. 7


Fig. 8


Fig. 9



Fig. 10

Eliminating air bubbles

If air bubbles should build up during application, they can be opened with a pointed instrument, e.g. a probe, and filled. You can also inject additional Orthocryl® LC from underneath the air bubble so that it moves towards the surface and disintegrates.

Application of colored Orthocryl® LC

In order to guarantee secure hardening of the color layers, we recommend applying a layer of clear Orthocryl® LC first.

Orthocryl® Disco glitter

Spray a small amount of Orthocryl® Disco glitter on the first layer of Orthocryl® LC directly from the spray bottle. Then, cover it with another layer of Orthocryl® LC.

Application of decals

Place the cut decal on the first Orthocryl® LC layer, at about the height of the wire retentions. Then, cover it with another layer of Orthocryl® LC.

5.1. Application of the Orthocryl® LC bite paste

Example: activator

After the upper and lower model have been placed in the fixator by means of an activating bite*, both parts of the fixator are separated and models are prepared as described in section 5.

Orthocryl® LC is then applied to the palatal side of the upper jaw and to the lingual side of the lower jaw. The Orthocryl® LC bite paste is then taken from the jar with a spatula or a wax knife and applied onto the occlusal or the incisal part of the lower model. To ensure the material is entering all areas of the structure (bite ramp) and to avoid air bubbles, make sure you apply an excess amount of bite paste.

Following this, both parts of the fixator are put together again, the excess material is removed, if necessary, and the closed fixator is placed in the light-curing unit.


*There is enough space in the Solidilite V unit from Shofu for average-sized models if they are placed in the fixator 072-004-00 in a space-saving manner (close to telescopic bars).

6. Polymerization

The specifications on polymerization times refer to the recommended light-curing unit Solidilite V from SHOFU Dental GmbH, Germany. Please observe the instructions for use for this unit.

The rotary plate in the curing chamber is height-adjustable. During polymerization of the described expansion plate, the rotary plate should be at the highest position possible according to the instructions for use of the light-curing unit Solidilite V.

The polymerization of Orthocryl® LC in the light-curing unit takes place at a wave length of 400 – 550 nm and a light output of 4 x 150 watt (600 watt). The time of polymerization for a single jaw appliance is 3 minutes and 6 minutes for bimaxillary appliances (Fig. 9). Afterwards, carefully remove the appliance from the model (Fig. 10) and cure the appliance for further 1 – 3 minutes with the basal side up.

Caution!

Remove the air-inhibition layer (smear layer) on the polymerized appliance, see point 7.

7. Removing inhibition layer

Completely remove the air-inhibition layer (smear layer). Apply Orthocryl® LC cleaning liquid (REF 160-430-00) on a piece of fabric and remove the smear layer.

Caution!

To eliminate any kind of risk for the patient, it is essential that the air-inhibition layer (smear layer) is completely removed (see safety data sheet). Direct application in the mouth is not permitted!

8. Finishing

When finishing Orthocryl® LC appliances, proceed as with Orthocryl® cold polymer appliances. However, do not use cross-cut hard metal burs. Hard metal burs with transversal toothings or burs for soft acrylics are better suited.

Pre-polish the appliance with silicone polishers or sand paper, and then with pumice powder and polishing brushes. A high polish finish can be achieved with an Edelweiss acrylic polishing material (REF 190-100-00 or REF 190-200-00) and a cloth polishing brush.

Note:

Due to stabilizing fillers, Orthocryl® LC appliances exhibit a slight milky haze.

9. Composition

See safety data sheet.

10. Safety instructions

For curing Orthocryl® LC, use light-curing units that emit light at a wavelength of 400 – 550 nm and have a light output of 4 x 150 watt (600 watt).

In order to avoid direct skin contact with Orthocryl® LC during processing, we recommend wearing protective gloves (e.g. Dermatril® disposable nitrile gloves, REF 162-907-00, REF 162-908-00, REF 162-909-00), see safety data sheet.

Direct application of Orthocryl® LC in the mouth is not permitted!

11. Storage and shelf life

Since Orthocryl® LC is a light-curing acrylic, close cartridges after use and store away from light.

The shelf life in the closed original package is 30 months from the date of fabrication:

- Protect from contamination and sunlight as they may cause premature polymerization.

12. Cleaning and disinfecting Orthocryl® LC appliances

Appliances should be cleaned regularly and stored in water when not in use to prevent problems with hygiene. This procedure prevents crack corrosion or brown discoloration of the wire components.

For professional cleaning in practices, we recommend using the activeblue high performance cleaning unit (REF 098-000-00). In combination with activefluid cleaning liquid (REF 150-020-00), specifically developed for orthodontic appliances, fast and effective cleaning of plaque and tartar is ensured, especially in areas that cannot be cleaned with a normal tooth brush.

13. Additional tips

Should the user and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the user and/or the patient is resident is informed accordingly.

14. Notes on disposal

When disposing of the product, please adhere to the national regulations that are valid for your region and - if available - the details outlined in the safety data sheets.

**Please refer to the label. Additional information can be found at www.dentaurum.com
(Explanation of the symbols REF 989-313-00).**

15. Delivery program

Orthocryl® LC starter set

REF 160-400-00

Content: Orthocryl® LC, 2 x 30 g of each color (clear, transparent pink, blue), 1 injector, 50 cartridge tips and step-by-step instructions

Orthocryl® LC starter set, clear

REF 160-400-20

Content: Orthocryl® LC, clear, 1 injector, 50 cartridge tips and step-by-step instructions

Orthocryl® LC single cartridges and accessories

Designation	Quantity		Included in:	
			starter set	starter set, clear
Orthocryl® LC, clear	30 g	REF 160-401-00	2 x	6 x
Orthocryl® LC, transparent pink	30 g	REF 160-402-00	2 x	
Orthocryl® LC, blue	30 g	REF 160-404-00	2 x	
Bite paste for Orthocryl® LC	50 g	REF 160-410-00		
Injector for Orthocryl® LC	1	REF 160-420-00	•	•
Cartridge tips for Orthocryl® LC	50	REF 160-421-00	•	•
Cleaning liquid for Orthocryl® LC	100 ml	REF 160-430-00		
Separating medium	1000 ml	REF 162-800-00		
3D separating medium	100 ml	REF 162-801-00		

16. Quality

Dentaurum ensures faultless quality of its products. These recommendations are based upon Dentaurum's own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as we, Dentaurum, have no influence on the processing on site.

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