Fig. 1 Fig. 2 Fig. 3

Fig. 4 Fig. 5

Fig. 6


Fig. 7


Fig. 8


Fig. 9 Fig. 10


Fig. 11 Fig. 12


Fig. 13 Fig. 14


Fig. 15 Fig. 16


Fig. 17 Fig. 18


Fig. 19

Fig. 20

ES

Estimado cliente:

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Dear customer,

Thank you for choosing a quality product from Denta... Please read these Instructions for use to ensure you use this product in a safe and simple way and that you and your patients gain full benefit.

1. Manufacturer

Denta... GmbH & Co. KG, Turnstr. 31 | 75228 Ispringen | Germany

2. General device description

The flexible forward thrust double plate system acc. to Schaneng differs from other systems by means of the flexibility of the intermaxillary guide elements. Each of the guiding ridges, which are connected to one another, is made of special 1.0 mm round spring wire.

3. Intended purpose

The elastic forward thrust double plate system acc. to Schaneng (removable orthodontic appliance) is used to correct the maxillo-mandibular relationship in patients with a skeletal malocclusion class III.

4. Intended users

All products are intended for dental professionals only.

5. Indications

Correction of the maxillo-mandibular relationship, in particular for a class II malocclusion. Transverse extension of the upper jaw or the lower jaw is also possible.

6. Contraindications

Stainless steel expansion screws must not be used in patients with a nickel and/or chromium allergy. Warnings are found on the respective product labels.

7. Patient target group

All patients for whom a removable appliance with a forward thrust double plate appliance can be used as part of orthodontic treatment.

8. Fabrication of the flexible forward thrust double plate

The first stage in making this flexible appliance begins by joining the two models with the construction bite in the fixator (Fig. 6). Elements for retention and movement, along with diverse screws (except for the upper jaw screw with guide ridges), can be attached to the models in the same manner...

The mounting template and the plastic adapter are required to set the angle of the inclined plane at 60° on the mandibular plate. The mounting template is laid on the lower jaw model parallel to the occlusal plane. At the side, the 60° guide surfaces of the template are positioned adjacent to the screw holder and tightly against the screw body.

The plastic adapter is now pushed over the guide surfaces of the template and attached to the screw holder by means of sticky wax (Fig. 10). When the mounting template has been removed, the position and distance of the plastic adapter in relation to the model and the screw can be seen and checked (Fig. 11). The mandibular plate is made in the usual way with the salt-and-pepper technique, polymerized in the pressure vessel and the plastic adapter is then removed (Fig. 12).

Before fitting the guide ridges of the maxillary plate onto the inclined plane, the mandibular plate must first be fully finished and polished, and the anterior region thinly covered in wax to avoid the possibility of the plates bonding together during polymerization (Fig. 13). Note: The loops on the guide ridges of the upper jaw screw are supplied waxed out. This ensures free movement/flexibility of the guide ridges in the acrylic. The wax layer should be checked and adjusted if necessary before applying acrylic (Fig. 14).

The guide ridges are now waxed into position parallel to the inclined plane. The screw itself should be placed as close to the palate as possible (Fig. 16). Following this, the acrylic is applied in the upper jaw area with the fixator closed (Fig. 17). The next stage is to finish the maxillary plate and to polish both plates. Ensure that the loops are completely free of acrylic and that the guide ridges are not damaged by the grinding tools (Fig. 18).

Important

The construction bite should always be taken on the patient taking the following into consideration: sagittal: max. 5 mm forward thrust, vertical: 3 mm opening in anterior region

9. Composition

The material composition is included in the orthodontics materials list; please refer to the orthodontics catalog or www.dentaurum.com.

10. Safety information

The product should not be used if there is a known allergic reaction to one or more of the material components. It is the dental professional's responsibility to decide on combining different products with regard to the compatibility of chemical, physical and mechanical properties.

11. Information for single use products

The product is intended for single use only. Reconditioning (recycling) of the product and its reuse is not permitted, as perfect functionality of the product cannot be guaranteed.

12. Further information

The patient brochure „Oral hygiene during orthodontic treatment“ (REF 989-328-20), available at www.dentaurum.com, gives information on how to keep this and similar appliances clean and maintained.

13. Product range

For the complete product range for the flexible forward thrust double plate system acc. to Schaneng, please see our orthodontics catalog or www.dentaurum.com

14. Quality information

Denta... ensures faultless quality of its products. These recommendations are based upon Denta...'s own experiences. The dental professional is solely responsible for the processing of the products.

15. Explanation of symbols used on the label

Please refer to the label. These Instructions for use and additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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