rema<sup>®</sup>  
TT **R**

Gebrauchsanweisung • Instructions for use •  
Mode d'emploi • Modo de empleo • Modalità d'uso

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## Dear Customer

Thank you for choosing a quality product from Dentaaurum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current instructions for use supplied with the product and stored in the internet at [www.dentaaurum.com](http://www.dentaaurum.com), even though you frequently use the same product.

## 1. Manufacturer

Dentaaurum GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

### Application of the rema® TT test specimens

The test specimen is used to determine the ideal concentration of rema® TT mixing liquid for the production of secondary telescope crowns.

The outer surface of the test specimen has been milled at an optimum angle of 0°. Sometimes laboratory milling machines show deviations from this ideal value (irregular motor, uneven bur). If this is the case for you, please consider the deviation in the actual end result!

#### Tip

The test specimen can be individualised if desired, i.e.

- by integrating a notch to act as a rotational stop.
- by adding an acrylic or plaster base.
- by re-milling with your own machine to check if your milling machine/bur actually produces a milled angle of 0°.

#### Application

Fabricate a secondary crown in modelling acrylic (e.g. Pattern Resin\*\*), vacuum-formed foil or wax.




Cover the surface with a thin layer of wax.

Invest according to the rema® TT instructions for use.

#### Fit

Too tight ➡ increase the mixing liquid concentration.

Too loose ➡ reduce the mixing liquid concentration.

\* registered trademark of the GC corporation

## 2. General description

Phosphate bonded speed investment material for the double crown technique using non-precious metal alloys.

The special advantages of this investment material are a high degree of material expansion, the freedom of choice concerning pre-heating methods, without having to forfeit the fitting accuracy, and consistently high surface smoothness. In order to benefit from the high quality casting results, it is necessary to observe the following instructions for use carefully.

## 3. Safety information

When using the speed heating method, the furnace door must remain closed for at least 15 minutes after the casting ring has been inserted. (Danger of burning due to spurting flame).

Investment materials contain quartz. Do not inhale the dust! Risk of lung damage (silicosis/lung cancer). Recommendation: Wear protective mask type FFP 2 – EN 149:2001. Cut open the bag using a pair of scissors and avoid generating dust when emptying into the mixing beaker. Fill the empty bag with water before scrunching-up.

Remove dust at the work place using a damp cloth.

In order to avoid dust generation during devestment, soak the ring in water until it is completely moist.

When sandblasting, use an extractor with a fine filter dust bag.

## 4. Delivery

<b>rema® TT Powder</b>	4 kg (25 x 160 g)	<b>REF 106-740-00</b>
<b>rema® TT Mixing Liquid</b>	1000 ml	<b>REF 106-707-00</b>
<b>rema® TT Test Specimen</b>	1 piece	<b>REF 319-740-00</b>

## 5. Service-Life

Powder – 36 months with dry storage

Liquid – 18 months (frost sensitive)

## 6. Material Properties

DIN EN ISO 15912 | Type **1**, class **1 2**

Setting begin	6 min*
Compressive strength	9 MPa*
Flowing characteristics	165 mm*
Thermic expansion	1.0 %*

\* Liquid 100 %

# 7. Application

**Storage**

**Wax wetting agent**

Only to be used on wax!  
Do not use on modelling acrylic!

**Casting ring system**

With the speed heating method,  
only use metal casting rings!

**Mixing ratio**

**160 g : 33 ml**

**Liquid concentration**

**50 % – 90 %**

**Mixing**

**90 sec.** Mixing time under vacuum.

**Working time**

**4 – 5 min**


**Investing**

Vibrate gently.

**Conventional heating method**

### Additional Information

Observe the ideal working temperature of 20 – 22 °C / 68 – 72 °F!

Only cool the liquid in case of higher room temperatures!

**Caution!** The total working temperature of powder and liquid is important. For best results, store in a climatic cabinet!

Temperature influences the rate of expansion.

Protect the mixing liquid from frost!

**Caution with deliveries in winter!**

When using a wax wetting agent, please use Lubrofilm® or Lubrofilm® plus! Please dry the pattern thoroughly! Without wax wetting agent, mix thoroughly under vacuum!

Avoid alcohol residues, apply thinly and blow-dry.

For casting ring sizes 3 and 6 please use a dry 2 mm ring liner (REF 127-251-00).

When using silicone or rubber casting rings, these must be soft and elastic. Line slightly with Vaseline!

The mixing ratio must be observed accurately!

Use a clean measuring beaker!

**See separate mixing table!**

**When using the speed heating method, increase the mixing liquid concentration!**

First put liquid, then powder in the mixing beaker. Mix by hand thoroughly first. Use a clean mixing beaker! No plaster residues!

Mixing speed and mixing blade geometry influence the mixture and therefore the expansion!

The ideal mixing speed is approx. 360 min<sup>-1</sup>.

Measured at room temperature 20 – 22 °C / 68 – 72 °F.

Cool the liquid at higher temperatures.

Set the vibrator to a low frequency and do not vibrate for too long.


Place the ring liner so that it is flush with the top edge of the ring. Fill the ring completely up to the ring liner. Roughen the upper surface of the investment material after setting!

Setting time: 40 min  
Ideal heating speed: 5 °C / 41 °F/min  
Holding time at 250 °C/h  
End temperature: 900 – 950 °C / 1652 – 1742 °F  
Holding time at end temperature: 40 – 60 min according to ring size


The base of the furnace should be corrugated. Place the ring into the furnace with the opening at the bottom! Protect the ring from drying out over the weekend! (seal the ring with wax)

**Reduced speed heating**

**Speed heating**

**Cooling down / devesting** Allow to cool down slowly to room temperature.

## 8. Recommended dilution of the mixing liquid

Indication	Premolar	
	Speed heating (metal casting ring)	Conventional (ring-free)
Primary crown / wax	<b>70 %</b>	<b>60 %</b>
Secondary crown (Pattern; Folident)		
Telescope 0° *	<b>70 %</b>	<b>60 %</b>
Conical 2°/4° *	<b>60 %</b>	<b>50 %</b>
For a molar unit, reduce the concentration by approx. 5 %!		
* An optimal fit is achieved, when the secondary crown covers two-thirds of the primary crown after the cast.		

**Tip:** The above values represent the results attained in the Dentaurum test lab. In individual cases these results can be influenced by different working methods, modelling acrylics and temperatures.

**Ring-free investing!**

Setting time: 30 min

Place ring in a pre-heated furnace at 800 °C/1472 °F!

Residence time: 15 min

Heat up consistently to 900 – 950 °C/1652 – 1742 °F!

Holding time at end temperature: 60 min

**Only use metal casting rings!**

Once the ring has been placed in the furnace, do not open until at least 15 min have elapsed – danger of instant combustion!

Pull the ring liner so that it is flush with the upper edge of the ring. Roughen the surface of the upper investment side with a knife. After 30 min place the ring immediately into the furnace at the end temperature! (900 – 950 °C/1652 – 1742 °F)

Holding time: 60 min

Holding time for larger rings (6 – 9): 60 – 90 min

When using plastic parts, always cover with a layer of wax!

Moisten casting rings and carefully devest.

Sandblast with glass beads or aluminium oxide at 50 µm.

Do not use a hammer to devest!

## 9. Mixing table

Concentration of mixing liquid	160 g : 33 ml	
	Mixing liquid concentration (ml)	Distilled water
90 %	30	3
80 %	26	7
70 %	23	10
60 %	20	13
50 %	17	16

## 10. Quality information

Dentaurum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

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