# NACERA<sup>®</sup>

## PEARL Q<sup>3</sup> MULTI-SHADE

GEBRAUCHSHINWEISE DE  
INSTRUCTIONS FOR USE EN  
ISTRUZIONI PER L'USO IT  
INSTRUCTIONS D'UTILISATION FR  
INSTRUCCIONES DE USO ES  
ИНСТРУКЦИЯ ПО ИСПОЛЬЗОВАНИЮ RU  
取扱説明書 JP

## Nacera® Pearl Q<sup>3</sup> Multi-Shade

DE Gebrauchsanweisung	4
EN Instructions for use	8
IT Istruzioni per l'uso	12
FR Instructions d'utilisation	16
ES Instrucciones de uso	21
RU Инструкция по использованию	26
JP 取扱説明書	31

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## EN Instructions for use

### Safety notes

Please read the instructions for use carefully before removing the blank (white made of zirconium oxide) from the packaging. They contain important information which is essential to ensure both the correct processing and the safety of patients and users.

Blanks made of **Nacera® Pearl Q<sup>3</sup> Multi-Shade** are manufactured and tested according to highest quality standards. In order to guarantee this level of quality during subsequent processing, the procedures described in the following must absolutely be adhered to.

### Warning

When using this product, the dentist must consider possible interaction between this medical product and other products already in place in the patient's mouth.

### Danger

The processing of blanks and finally-sintered mounts creates dust which might damage the lungs as well as irritating the eyes and skin. Therefore, processing may only be done provided that the extractor system works correctly as well as wearing safety goggles and a filtering facepiece. Please forward all of the information above to the dentist if you process this product for use in special designs within the framework of Directive 93/42 EEC on Medical Devices.

Please also note the instructions for use and safety data sheets.

### General handling notes

Blanks made of **Nacera® Pearl Q<sup>3</sup> Multi-Shade** are supplied in a pre-sintered white-fired condition. They have limited stability or strength and a residual porosity. Therefore, these blanks must be handled with care.

The blanks must be stored in their original packaging (0 - 90% air humidity), at temperatures between 10 °C and 50 °C. They must not be subjected to impact or vibration. Contamination must absolutely be avoided. Make sure that the mounts are only handled with dry and clean hands or gloves and that they are under no circumstances contaminated with liquids (such as adhesives or marker pens).

If tooth stumps and metal mounts are heavily discoloured, **Nacera® Pearl Q<sup>3</sup> Multi-Shade** should not be used due to its high translucence. In this case, we recommend **Nacera® Pearl Multi-Shade**.

Coolants will reduce the material's translucency.

### Indication

**Nacera® Pearl Q<sup>3</sup> Multi-Shade** is zirconium oxide 6Y-PSZ (Partly Stabilized Zirconia) for the production of fixed dental prostheses. The material is suitable for single crowns and/or bridges in the front and side tooth areas consisting of up to 3 units.

In the posterior region, the span between the pillars must not exceed a premolar width. For bruxism, the vertical dimension must be observed. In these cases, please consult with the dentist.

### Contraindications

In the event that there is an insufficient vertical dimension and that the preparation is unsuitable for an all-ceramic restoration, an alternative material must be chosen. Inlay bridges, enossal implants and root posts are other contraindications.

### Preparation

Recommendations for the preparation include a marked chamfer or a rounded step.

- Minimum cutting depth at the preparation border 1 mm
- 1.5 to 2 mm of substance removed occlusal/incisal
- Edge radius 0.7 mm
- Preparation angle 6° to 8°

For bridge constructions, observe the parallelism and avoid negative steps. In general, please observe the notes in the specialist literature.

### Notes on construction

Please place the construction in the layers of the **Nacera® Pearl Q³ Multi-Shade** blanks in so that the required colour/colour intensity is achieved. For more detailed instructions, please order the colour guide **Nacera® Pearl Q³ Multi-Shade** or download it in the download area at [www.nacera.de](http://www.nacera.de).

**Wall thickness Nacera® Pearl Q³ Multi-Shade:** For all-ceramic dentures, the wall thickness of the crowns in the sintered state must not be less than the corresponding values of table\*. The minimum wall thickness must always be observed even after any adjustment of the mount by the dentist. The geometry of the connectors is of the greatest importance for the crack resistance of mounts made of zirconium oxide. Therefore, the cross-section of the connectors must as large as possible.

For static reasons, the height of the connector is especially important. In order to avoid so-called „chipping“, the mounts should correspond to a reduced, fully anatomic shape of the crowns and bridges to ensure that the applied ceramic has the best possible support.

According to manufacturer specifications, the veneering ceramic should be layered at a layer thickness between 0.7 and 1.5 mm, but < 2 mm. Furthermore, the mounts must be designed according to the general guidelines of digital dental technology.

### \*Wall and edge thickness

	Single constructions	Bridges
Wall thickness minimum	0.4 mm	0.5 mm
Edge thicknesses minimum	0.2 mm	0.2 mm

## Mount dimensions

Number of intermediate segments	1
Connector cross section	9 mm <sup>2</sup>

## Milling strategy

In your CAM software, for **Nacera® Pearl Q<sup>3</sup> Multi-Shade**, please select the "highly translucent strategy".

## Processing

Only those machines and tools which are approved for processing pre-sintered blanks made of zirconium oxide may be used for manufacturing dentures made of **Nacera® Pearl Q<sup>3</sup> Multi-Shade**.

Caution! For **Nacera® Pearl Q<sup>3</sup> Multi-Shade**, the printed side is the occlusal side.

If possible, neither coolant nor compressed air should be used during processing. Please observe the instructions for use of your milling machine as well as the parameters of the CAD/CAM software.

Please separate the finished milled parts from the blanks carefully using suitable tools. Then thin the edges which have been thickened by the machine and remove the tapping. You can now carry out a cut-back manually if this has not already been done in the CAD.

## Visual inspection

Before the milled mounts are processed further, they must

be inspected for the following faults:


- Shiny areas on the surface (indicating a worn milling cutter)
- Discolouration (see also the user instructions for **Nacera® Clean**)
- Material spalling (due to the milling strategy and milling cutter)
- Cracks

Faulty mounts must not be processed further.


## Final sintering

The final sintering, also called dense sintering, ensures that the milled workpiece (mount) is given its final properties. During the final sintering, the mount shrinks by a precisely defined factor. In order to do this, the described procedure must absolutely be adhered to. Sintering can take place in all common dental sintering furnaces which are approved for the sintering of mounts made of zirconium oxide. Since zirconium oxide is known to be a poor thermal conductor, it is recommended that you slowly heat the mounts to the required temperature (see sintering graph) and cool them down equally slowly. Sintering supports should be thinned out or milled as rods. The ideal translucence is achieved during a sintering procedure at 1450 °C.

To perform just one sintering procedure with other Nacera products, **Nacera® Pearl Q<sup>3</sup> Multi-Shade** can be sintered at 1500 °C with slightly reduced translucence (approx. 2 %).

*Sinter curve with ideal translucence.*

*Sinter curve with slightly reduced translucence of approx. 2%.*

### Further processing of the finally-sintered restorations

If possible, the finally-sintered units should not be processed mechanically. However, if this is unavoidable, only water-cooled diamond-fitted tools in correct working order may be used. Ideally, the interdental connection points should not be polished. In principle, the basal polishing of these connection points must always be avoided for stability reasons (predetermined breaking points). Also in the case of abutments, sharp edges should be avoided if possible. Instead, rounding is to be aimed for.

### Veneering

Veneering is done with commonly available veneer ceramics that are approved for zirconium oxide

according to manufacturer specifications.

### Individualisation

For individualising the restoration, painting techniques as well as cut-back and layering techniques, or a combination of both, are suitable.

### Grinding

To protect the antagonists (abrasion) and for material engineering reasons, ground contact points and surfaces must either be high-gloss polished after the try-in and/or glaze-fired. To restore the physiological characteristics of the surface, polishing with Nacera® Shine Zr is necessary.

### Incorporating

The inner surfaces of the restoration can be fitted with a mechanical retention by means of blasting with  $\text{Al}_2\text{O}_3$  50  $\mu\text{m}$  at max. 1.5 bar. Before the retention/fastening system is inserted, the blasted interior surface should be cleaned with alcohol. State-of-the-art dental technology prefers self-adhesive and adhesive fastening over conventional fastening.

The fastening logs of the cement manufacturers must be observed.

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DOCERAM Medical Ceramics GmbH | Hesslingsweg 65-67 | D-44309 Dortmund