## remanium® star

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### EN

## remanium® star

**CoCrW - Dental alloy based on cobalt for metal-ceramic bonding, type 5 acc. to DIN EN ISO 22674 / DIN EN ISO 9693 for melting and casting**

**Composition (% by mass)**

Co	Cr	W	Si
60.5	28.0	9.0	1.5

Additional elements < 1%: Mn, N, Nb

This product contains cobalt.

Free from nickel, beryllium, lead, cadmium and iron.

Yield strength*	R <sub>p0,2</sub>	620 MPa
Tensile strength*	R <sub>m</sub>	845 MPa
Hardness*	H	280 HV10
Elongation at rupture*	A <sub>5</sub>	10.2%
Modulus of elasticity*	E	190 GPa
Density	ρ	8.6 g/cm³
Solidus temperature	T <sub>s</sub>	1320 <span> </span> °C / 2408 <span> </span> °F
Liquidus temperature	T <sub>l</sub>	1420 <span> </span> °C / 2588 <span> </span> °F
Pouring temperature (recommended)	T <sub>g</sub>	1520 <span> </span> °C / 2768 <span> </span> °F
CTE (25 <span> </span> °C - 500 <span> </span> °C / 77 <span> </span> °F - 932 <span> </span> °F)	α	14.1 x 10 <sup>-6</sup> K <sup>-1</sup>
Weight per casting	m	approx. 6 g

\* Variations of ± 10 % are possible depending on the batch.

**Dear customer,**

Thank you for choosing a quality product from Dentaaurum. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaaurum.com, even though you may frequently use the same product.

**1. Manufacturer**

Dentaaurum GmbH & Co. KG I Turnstraße 31 I 75228 Ispringen I Germany

**2. General product description**

Non-precious alloy for the crown and bridge technique.

**3. Area of application**

remanium® star is a cobalt-based cast alloy. It is particularly suitable for the production of crowns and bridges, double crowns, implant-supported superstructures, adhesive bridges, secondary parts for cast partial dentures (combination work) and for metal ceramics.

**4. Wax-up**

Construction must be in line with dental regulations in order to fulfil mechanical requirements. Use only materials (e.g. wax, plastic) that burn without residue. Detailed information can be found in the processing instructions.

**5. Investing**

Investment materials for non-precious alloys with high overall expansion must be used, such as rema® CC, Castorit® all speed, Castorit®-super C, Trivest and rema® TT or rema® CAD/Vest (Dentaaurum).

**6. Melting and casting**

Only suitable ceramic-based melting crucibles should be used for remanium® alloys. Care must be taken to ensure that for each alloy separate melting crucibles are used which are not contaminated and/ or damaged.The use of new material is a prerequisite for ensuring clear batch traceability. Suitable melting and casting processes are high-frequency melting with vacuum pressure casting or centrifugal casting machines, arc melting with vacuum pressure casting or open melting. After the breakup of the melt, cast off immediately to avoid overheating of the alloy. Following the casting process, the casting ring is carefully removed and should cool slowly in the air.

**7. Finishing**

The finishing of the framework and the surfaces should be carried out by means of finely-toothed hard metal burs, ceramic-bonded grinding points or diamond grinding points.

**8. Welding**

Suitable welding processes such as laser welding and TIG impulse welding produce a solder-free, mechanically strong and corrosion-resistant joint. It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters of each individual laser machine. Suitable welding wires are available as filler material for all remanium® alloys, e.g. CoCr welding wire 0.25 mm (REF 528-215-10), 0.35 mm (REF 528-210-10) and 0.5 mm (REF 528-200-10).

**9. Soldering**

Soldering should be avoided if possible to reduce the number of materials. If, however, it is necessary to solder, please use a suitable solder and flux which are appropriate for the composition and melting temperature of the alloy used., e.g. CoCrMo-Sold 1 (REF 102-306-00) in conjunction with rema® Flux 1 (REF 102-304-00).

**10. Ceramic veneering**

remanium® bonding alloys can be veneered with all suitable ceramic masses, e.g. ceraMotion® Me (Dentaaurum), as long as the CTE is suitable for the alloy. The highest recommended firing temperature is 980 °C / 1796 °F. The ceramic should cool down over a long period in the furnace unless the ceramic manufacturer has other specifications.

**11. Contraindications and adverse reactions**

Signs of intolerance to non-precious remanium® alloys are extremely rare if the manufacturer’s Instructions for use are adhered to. If the patient has a proven allergy against any component within the alloy, this alloy must not be used for safety reasons. There have been individual reportings of local irritations which were electrochemically induced. If various alloys have been used, it is possible that galvanic effects may occur.

Exposure to dust and vapors may cause irritation to the eyes and/or respiratory tract. Individual alloy components may have a carcinogenic effect. Please adhere to warnings and safety data sheets. In the State of California elements of this alloy are known to be carcinogenic.

**12. Warnings and precautions**

Metal dust may be produced as the alloy is mechanically processed. The mechanical process will also create heat and possibly also burs. Processed material may therefore be hot and/or have sharp edges. Processing at temperatures above the solidus temperature may generate vapors. The workplace should be fitted with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and respiratory protection equipment are worn. Should there be contact with the skin, please wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, please rinse the eye with a gentle stream of water or saline solution for at least 15 minutes. Please consult a doctor if you are exposed to individual elements.

**Please observe the relevant safety data sheet (see www.dentaaurum.com).**

Interproximal or occlusal contact with other metallic elements and the use of different types of alloys in the same oral cavity should be avoided to prevent galvanic reactions or local irritations in individual cases. The alloy should not be ground or polished in the patient’s mouth. There is no scientific evidence on the safety or efficacy of treatment of pregnant women or nursing mothers or children.

**13. Reuse**

The alloy may be changed in its composition and properties if reused. Its reuse is therefore not permitted. This is the only way to ensure clear batch traceability.

**14. Conditions for storing**

Normal changes to the surrounding conditions (e.g. temperature, pressure, light) have no negative impact on the product characteristics of the remanium® alloys.

**15. Disposal**

Please adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

**16. Quality**

Dentaaurum ensures faultless quality of its products. These recommendations are based upon Dentaaurum’s own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as we, Dentaaurum, have no influence on the processing on site.

**17. Further information**

Should the user and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the user and/or the patient is resident is informed accordingly

The SSCP is available at https://ec.europa.eu/tools/eudamed.

**18. General information**

These Instructions for use contain the main steps for processing the non-precious remanium® alloys. For more information, see remanium®, Non-precious metal alloy for cast partial denture technique REF 989-448-00.

**19. Explanation of symbols used on the label**

⚠ Please refer to the label. Additional information can be found at www.dentaaurum.com (Explanation of symbols REF 989-313-00).

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