**Estimado cliente**

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| EN |
|----|
|----|

### Dear customer,

Thank you for choosing a quality product from Dentaurum. Please read these Instructions for use to ensure you use this product in a safe and simple way and that you and your patients gain full benefit. Should you have any questions or ideas, please contact your local representative.

### 1. Manufacturer

Dentaurum GmbH & Co. KG
Turnstr. 31 | 75228 Ispringen | Germany

### 2. General product description

To use on the upper or lower jaw archwire of a multi-bracket appliance. The Herbst II retention hinge is fitted with a special connecting element with which the hinges can be attached to the archwires already fitted on the multi-bracket appliance. This enables combined use of the Herbst hinge and a multi-bracket appliance without patient compliance.

### 3. Intended purpose

Correction of distoclusions.

### 4. Intended users

All products are intended for dental professionals only.

### 5. Indications

Distocclusion/class II (unilateral/bilateral), molar distalization, gap closure (aplasia), temporomandibular dysfunction, dentoalveolar compensation of occlusion (elastics substitution).

### 6. Contraindications

The Herbst appliance made of stainless steel must not be used for patients with a nickel allergy. Warnings are found on the respective product labels. Also see 9. Composition.

### 7. Patient target group

All patients for whom a fixed appliance can be placed as part of ortho-dontic treatment.

### 8. Processing

To fit the Herbst II hinge, the multi-bracket appliance must already be fitted with the strongest steel archwire possible for this treatment technique.

The connecting element (nut) of the hinge has a slot via which it is placed on the archwire on the multi-bracket appliance and fixed with the supplied screw. The hinges are connected this way to the upper jaw, first molar mesially and to the lower jaw on the archwire between the cuspid and first bicuspid.

Intra-oral procedure:

After the nut has been placed on the steel archwire, the distance is measu-red between the upper and the lower nut and the hinges are shortened accordingly. Then the hinges are fixed onto the nuts with the screws. The two longer screws are usually used in the upper jaw and the shorter ones in the lower jaw. Depending on the situation, the screws can be used the other way round, for example if the hinge is in the way of the brackets. If the archwire of the multi-bracket appliance is to be changed, the nuts must first be unscrewed. The hinge is then refitted as described.

If required, the Herbst II appliance can be activated with various spacer rings (see product range).

### 9. Composition

The Herbst II hinge is made of the following materials: DIN 1.4305 and DIN 1.4542. The material composition is included in the orthodontics materials list; please refer to the catalog or www.dentaurum.com. The SSCP is available at https://ec.europa.eu/tools/eudamed and www.dentaurum.com.

### 10. Safety information

The fit and the tightness of the screws should be regularly checked. If necessary, tighten the screws and secure them with a special screw securing adhesive, e. g. NOBIL GRIP.

The loops of the Herbst II retention hinge must not be bent. It is not intended that the loops can be ground.

The product should not be used if there is a known allergic reaction to one or more of the material components. If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and if necessary must be disposed of.

### 11. Information for single use products

The product is intended for single use only. Reconditioning (recycling) of the product and its reuse is not permitted, as perfect functionality of the product cannot be guaranteed.

### 12. Further information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the dental professional and/ or the patient is resident are informed accordingly.

### 13. Notes on disposal

When disposing of the product, please adhere to the national regulati-ons that are valid for your region and - if available - the details outlined in the safety data sheets.

### 14. Product range

For the complete range of Herbst products, please see the orthodontic catalog or visit www.dentaurum.com.

### 15. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as we,Dentaurum, have no influence on the processing on site.

### 16. Explanation of symbols used on the label

 **Please refer to the label. The Instructions for use and additional information can be found at [www.dentaurum.com](http://www.dentaurum.com) (Explanation of symbols REF 989-313-00).**

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