DE

remanium® GM 900

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EN

remanium® GM 900

CoCrMo – Dental alloy based on cobalt for cast partial denture technique, type 5 acc. to DIN EN ISO 22674 for melting and casting

Co	Cr	Mo	Si	Ta
62.0	30.0	5.5	1.2	1.2

Additional elements < 1%: N

This product contains cobalt.

Free from nickel, beryllium, lead, cadmium and iron.

Technical data

Yield strength*	R _{0,2}	700 MPa
Tensile strength*	R _m	910 MPa
Hardness*	H	340 HV10
Elongation at rupture*	A ₅	5%
Modulus of elasticity*	E	210 GPa
Density	ρ	8.2 g/cm³
Solidus temperature	T _S	1360 °C / 2480 °F
Liquidus temperature	T _L	1410 °C / 2570 °F
Casting temperature (recommended)	T _G	1510 °C / 2750 °F
Weight per casting	m	approx. 6 g

* Variations of ± 10 % are possible depending on the batch.

Dear customer,

Thank you for choosing a quality product from Dentaurum. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG · Turnstraße 31 · 75228 Ispringen · Germany

2. General product description

Non-precious alloy for cast partial denture technique.

3. Area of application

remanium® GM 900 is a cobalt-based cast alloy. It is particularly suitable for the manufacture of removable restorations with a rigid construction and a particularly high strength in the area of cast partial denture and combi technique.

4. Wax-up

Construction must be in line with dental regulations in order to fulfil mechanical requirements. Use only materials (e.g. wax, plastic) that burn without residue. Detailed information can be found in the processing instructions.

5. Investing

Investment materials for non-precious alloys with optimally adapted expansion must be used, such as rema® Exakt, rema® Exakt F and rema® dynamic S or rema® CAD/Vest (Dentaurum).

6. Melting and casting

Only suitable ceramic-based melting crucibles should be used for remanium® alloys. Care must be taken to ensure that for each alloy separate melting crucibles are used which are not contaminated and/ or damaged. The use of new material is a prerequisite for ensuring clear batch traceability. Suitable melting and casting processes are high-frequency melting with vacuum pressure casting or centrifugal casting machines, arc melting with vacuum pressure casting or open melting. After the breakup of the melt, cast off immediately to avoid overheating of the alloy. Following the casting process, the casting ring is carefully removed and should cool slowly in the air.

7. Finishing

The finishing of the framework and the surfaces should be carried out by means of finely-toothed hard metal burs, ceramic-bonded grinding points or diamond grinding points.

8. Welding

Suitable welding processes such as laser welding and TIG impulse welding produce a solder-free, mechanically strong and corrosion-resistant joint. It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters of each individual laser machine. Suitable welding wires are available as filler material for all remanium® alloys, e.g. CoCr welding wire 0.25 mm (REF 528-215-10), 0.35 mm (REF 528-210-10) and 0.5 mm (REF 528-200-10).

9. Soldering

Soldering should be avoided if possible to reduce the number of materials. If, however, it is necessary to solder, please use a suitable solder and flux which are appropriate for the composition and melting temperature of the alloy used., e.g. CoCrMo-Sold 1 (REF 102-306-00) in conjunction with rema® Flux 1 (REF 102-304-00).

10. Contraindications and side effects

Signs of intolerance to non-precious remanium® alloys are extremely rare if the manufacturer’s Instructions for use are adhered to. If the patient has a proven allergy against any component within the alloy, this alloy must not be used for safety reasons. There have been individual reportings of local irritations which were electrochemically induced. If various alloys have been used, it is possible that galvanic effects may occur.

Exposure to dust and vapors may cause irritation of the eyes and/or respiratory tract. Individual alloy components may have a carcinogenic effect. Please adhere to warnings and safety data sheets. In the State of California elements of this alloy are known to be carcinogenic.

11. Warnings and precautions

Metal dust may be produced as the alloy is mechanically processed. The mechanical process will also create heat and possibly also burrs. Processed material may therefore be hot and/or have sharp edges. Processing at temperatures above the solidus temperature may generate vapors. The workplace should be fitted with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and respiratory protection equipment are worn. Should there be contact with the skin, please wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, please rinse the eye with a gentle stream of water or saline solution for at least 15 minutes. Please consult a doctor if you are exposed to individual elements.

Please observe the relevant safety data sheet (see www.dentaurum.com).

Interproximal or occlusal contact with other metallic elements and the use of different types of alloys in the same oral cavity should be avoided to prevent galvanic reactions or local irritations in individual cases. The alloy should not be ground or polished in the patient’s mouth. There is no scientific evidence on the safety or efficacy of treatment of pregnant women or nursing mothers or children.

12. Reuse

The alloy may be changed in its composition and properties if reused. Its reuse is therefore not permitted. This is the only way to ensure clear batch traceability.

13. Conditions for storing

Normal changes to the surrounding conditions (e.g. temperature, pressure, light) have no negative impact on the product characteristics of the remanium® alloys.

14. Disposal

Please adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

15. Quality

Dentaurum ensures faultless quality of its products. These recommendations are based upon Dentaurum’s own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as we, Dentaurum, have no influence on the processing on site.

16. General information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident are informed accordingly. The SSCP is available at https://ec.europa.eu/tools/eudamed and www.dentaurum.com. These Instructions for use contain the main steps for processing the non-precious remanium® alloys. For more information, see remanium®, Non-precious metal alloys for the cast partial denture technique REF 989-451-00.

For more information on our products and services, please visit www.dentaurum.com. You will find these Instructions for use and information on the current status at www.dentaurum.com/ifu.

17. Explanation of symbols used on the label

Please refer to the label. Additional information can be found at: www.dentaurum.com/ifu (Explanation of symbols REF 989-313-00).

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