Gebrauchsanweisung | Instructions for use |  
Mode d'emploi | Modo de empleo | Modalità d'uso

**Orthocryl<sup>®</sup>**

**D**  
DENTAURUM

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## Dear customer,

Thank you for choosing a quality product from Dentaureum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at [www.dentaureum.com](http://www.dentaureum.com), even though you may frequently use the same product.

## 1. Manufacturer

Dentaureum GmbH & Co. KG

Turnstr. 31 | 75228 Ispringen | Germany

## 2. General product description, area of application

Orthocryl®/Orthocryl® EQ are cold-curing resins based on PMMA, used for the fabrication of fixed and removable dental appliances, occlusal splints, anti-snoring devices etc.

All Orthocryl® powders (polymers), fluids (monomers) and color concentrates are approved as medical products Class IIa for application for fixed orthodontic appliances.

Accessories such as Orthocryl® Disco glitter and decals are intended to be used for removable orthodontic appliances only.

**Orthocryl®** attains optimal stability when used with the salt-and-pepper technique. It is also suitable for the doughing technique.

When used with the doughing technique, **Orthocryl® EQ** attains optimal flowability. It is also suitable for the salt-and-pepper technique.

### 3. Processing

It is mandatory to wear face protection/safety glasses (REF 162-920-00) and special protective gloves, e.g. Dermatril® disposable gloves made of special nitrile (REF 162-907-00, 162-908-00 and 162-909-00) to avoid direct skin contact when processing Orthocryl® components. Please refer also to the corresponding safety data sheets.

#### 3.1. Plaster model/plastic model

Plaster models must be vented during the first step in processing. Otherwise, air bubbles from the plaster model could lift the acrylic or penetrate into the acrylic. Models should therefore be placed in lukewarm water for 10–25 minutes, depending on the quality of the plaster, until no more air bubbles rise from the plaster.

This step is not necessary in the case of plastic models.

When working with plaster models, we recommend using a special separating medium (REF 162-800-00) in addition after soaking in order to attain a smooth, silky basal surface on the appliance and to avoid discolorations and plaster inclusions.

Plastic models **must** be given a layer of separation medium before Orthocryl® is applied in order to prevent the plastics from sticking. We recommend the 3D separating medium (REF 162-801-00). This is, however, not available in the USA.

##### 3.1.1 Application of the separating medium (REF 162-800-00) on plaster models

- Blow the soaked model carefully with compressed air.
- Apply the separating medium from the dosing bottle (REF 162-100-00) directly on the soaked, slightly damp model and spread evenly using a brush.

- Stand the model on its back so that excess fluid can drain away.
- After approx. 30 seconds, carefully blow away separating medium residues or droplets with compressed air and, if necessary, remove residues from screws and wire retentions with a dry brush. Application of the acrylic can then start immediately using the salt-and-pepper technique (3.3) or the doughing technique (3.4).

### **3.1.2 Application of the 3D separating medium (REF 162-801-00) on plastic models**

To fabricate dental models, we recommend using high quality plastic models that were made using DLP (Digital Light Processing). e. g. prime4me® base model (REF 430-110-10). Besides being highly accurate and therefore ensuring very precise fitting, these models allow for the usual process in the laboratory.

Due to the large variety of materials for model printing, it is essential to test the suitability of the 3D separating medium before its use.

- Apply the 3D separating medium evenly with a brush over the entire model before the expansion screw (or screws) and wire elements are finally fixed in place and allow the separating medium to dry thoroughly. To this end, apply the separating medium from the dosing bottle (REF 162-100-00) directly on the model and spread evenly using a brush.
- Stand the model on its back so that excess fluid can drain away.
- Allow the layer of separating medium to dry thoroughly.
- The separating medium is completely dry when it no longer feels sticky.
- Repeat the process in order to achieve an even separating layer everywhere.
- In order to accelerate the drying of the separating medium or to shorten the drying phases, careful blowing with compressed air or a household hair dryer can be helpful. For faster drying of the separating medium, the model can also be placed in a light curing unit for a few minutes\*.

- \* The time required in the light curing unit may vary depending on the size of the workpiece and the amount of separating medium. The drying time also depends on the type of device – with the Solidilite V unit from Shofu, up to 3 minutes should be enough.
- If required, block out undercuts in the tooth/jaw area and fix the expansion screw(s) and wire elements as usual on the dried plastic model.
  - Application of the acrylic can then start immediately using the salt-and-pepper technique (3.3) or the doughing technique (3.4).

### 3.2. Salt-and-pepper technique

Apply first a layer of powder (polymer), then a layer of fluid (monomer) and continue applying alternate layers of polymer and monomer in slight pendulum movements

Only apply as much fluid as the powder can soak up. The material should not flow away.

To avoid the formation of air pockets under the screws, apply a slightly larger amount of fluid to begin with so that the acrylic spreads under the screws.

Finish application with a layer of powder. Apply enough powder so that the top layer is dry. This ensures a precise fit and minimal shrink.

### 3.3. Doughing technique

Mix the Orthocryl®/Orthocryl® EQ in a silicone mixing bowl at a ratio of 2.5 : 1 (powder : fluid).

Allow to swell – Orthocryl®, 5–7 minutes, Orthocryl® EQ 3–5 minutes (possibly shorter if room temperature is high), then apply the acrylic dough directly on the model.

### **3.4. Processing Orthocryl® color concentrates**

Orthocryl® fluids can be colored individually using the Orthocryl® color concentrates. We recommend a mixing ratio of 1 : 10 to 1 : 50 (color concentrate : Orthocryl® fluid clear). Precise dosing can be achieved with the dispensing pipette (REF 162-101-00).

Our recommendations for standard colors:

red 1 : 10                      blue 1 : 50

### **3.5. Processing Orthocryl® Disco glitter**

Orthocryl® Disco glitter can be mixed with Orthocryl® powder to the intensity required. The amount of glitter used has no effect on the properties or the processing of Orthocryl®. However, we recommend adding only a small amount of glitter, for example a small heap on the end of a wax knife or a Le Cron instrument. Mixing ratio: 0.5 : 100 (glitter : powder).

### **3.6. Processing decals**

Decals have a particularly striking effect when used with transparent acrylic. The decals and the acrylic do not bind together chemically so it is important to embed them completely in the acrylic.

### **3.7. Processing Orthocryl® black & white**

The opaque coloring of Orthocryl® black & white should be taken into consideration during processing. Orthocryl® Disco glitter and decals lose their impact in opaque acrylic. In particular during the finishing phase, avoid grinding the screws and wire elements. The opacity of Orthocryl® black & white can be reduced by adding clear Orthocryl®.

Mixing 50 % clear polymer with 50 % black or white polymer will reduce the opacity so that wire elements, screws etc. become visible in the Orthocryl®, yet the acrylic will still appear black or white. You can also choose other mixing ratios to mix the polymers and to obtain more or less opacity.

## 4. Time required for hardening / polymerization

Pressure vessel: 2.2 bar

Temperature: 40–46 °C / 104–114.8 °F

Time: **not less than 15 minutes** (Orthocryl® EQ)

**not less than 20 minutes** (Orthocryl®)

Do not interrupt polymerization in the pressure vessel during the first 8 minutes (Orthocryl® EQ) or 10 minutes (Orthocryl®) to prevent the formation of micro bubbles caused by boiling.

**Instructions for use and safety information for the pressure vessel should be adhered to strictly.**

## 5. Finishing and polishing

Orthocryl® acrylics can be finished and polished with the usual mills and polishers. Protective equipment (face and eye protection) should be worn.

## 6. Avoiding residual monomer and the biocompatibility of Orthocryl® acrylics

If Orthocryl® acrylics are processed correctly, the remaining monomer content is 3% by weight. To reduce the amount of residual monomer, place the finished appliance – if possible – in water for three days before being used by the patient. This will reduce the amount of residual monomers to similar values as for heat-curing acrylic. If this step is not possible, advise the patient to place the appliance in water whenever it is not being worn.

Orthocryl® acrylics have undergone scientific testing. Information on scientific testing and a detailed brochure of the entire Orthocryl® product range with many practical tips on processing is available on request.

- „Contribution to the Biological Assessment of Orthodontic Acrylic Materials“ in Fortschritte der Kieferorthopädie 200, 61, 246-257 (No. 4).

## 7. Information on cleaning and disinfection in the laboratory

Orthocryl® acrylics can be cleaned with a tooth brush and tooth paste. Please refer also to "Oral hygiene during orthodontic treatment", REF 989-328-20.

Orthocryl® acrylics can be disinfected with standard disinfecting agents for impression materials (e.g. Impresept® / 3M Espe, MD 520 / Dürr Dental, Dentavon® / Schülke+Mayr or Mucalgin® / Merz Hygiene). Follow the manufacturer's instructions.

## 8. Composition

Monomer = Methyl methacrylate "MMA" (fluid)

Polymer = Polymethylmethacrylate "PMMA" (powder)

Please refer also to the safety data sheets which are delivered automatically to customers with every first order and which are available at any time on request.

## 9. Information on dangers and safety

Contains methyl methacrylate and ethylene glycol dimethacrylate.

Causes irritation to the skin.

Can cause allergic reactions.

Causes serious irritation to the eyes.

Can cause irritation to the respiratory tract.

Keep away from heat/sparks/open flames/hot surfaces. Do not smoke.

Wear protective gloves/protective clothing/eye protection/face protection.

Please refer also to the safety data sheets which are delivered automatically to customers with every first order and which are available at any time on request.

The product should not be used if there is a known allergic reaction to one or more of the material components.

## 10. Storage and shelf life

Powder (polymer): 10 years in closed container.

Fluid (monomer): 36 months from date of manufacture in originally sealed bottle.

Protect from contamination and sunlight as they may cause premature polymerization.

## 11. Disposal

Please take note of the special instructions for the disposal of non-hardened Orthocryl® components as detailed in the safety data sheets.

## 12. Information for patients on cleaning and disinfecting Orthocryl® appliances

Appliances should be cleaned regularly and stored in fresh water when not in use to prevent problems with hygiene. This procedure prevents crack corrosion or brown discoloration of the wire components.

For more information, please refer to the pamphlet "Oral hygiene during orthodontic treatment", REF 989-328-20. This can be downloaded at [www.dentaurum.com](http://www.dentaurum.com).

## 13. Additional tips

Should the user and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the user and/or the patient is resident is informed accordingly.

## 14. Notes on disposal

When disposing of the product, please adhere to the national regulations that are valid for your region and - if available - the details outlined in the safety data sheets.

## 15. Delivery program

Please refer to our Orthodontics Catalog and our Prosthetics Catalog for information on the products available in the Orthocryl® range.

## 16. Quality

Dentaurum ensures faultless quality of its products. These recommendations are based upon Dentaurum's own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as Dentaurum has no influence on the processing on site.

**Please refer to the label. Additional information can be found at [www.dentaurum.com](http://www.dentaurum.com) (Explanation of symbols REF 989-313-00).**

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