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Abbreviated instructions for use

Dear customer,

Thank you for choosing a quality product from Dentaurem.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurem.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurem GmbH & Co. KG · Turnstr. 31 · 75228 Ispringen · Germany

2. General device description

The tiologic® TWINFIT implant system contains specially coordinated instruments, abutments, and accessory components for the prosthetic restoration of the tiologic® TWINFIT implants. Only original components of the tiologic® TWINFIT implant system should be combined in accordance with the Instructions for use.

The abutment is selected by the dentist/dental technician to suit the individual patient’s needs.

The diameter of the implant and the gingival thickness determine the choice of abutment.

Selection aids with conical and platform connector geometry for the abutment series S, M and L are available to help choose the right abutment on the model (gingival height and angulation).

3. Area of application

The tiologic® TWINFIT CAD/CAM components are used for the transfer and manufacture of digital prosthetic solutions such as single-tooth and bridge restorations on tiologic® TWINFIT implants.

4. Medical indications/contraindications

CAD/CAM components tiologic® TWINFIT

Indications

- Individual crown and bridgework (cemented, individually screw-retained)
- Abutment on implant ø 3.3 mm only in the area of the upper lateral incisors and lower incisors

Contraindications.

- Primary bracing of abutments

- Correction of angulation exceeding 20°, or exceeding 15° in the case of implants ø 3.3 mm

- Cast-on technique

5. Information on use, processing, fitting and insertion

The tiologic® TWINFIT CAD/CAM components are designed for the transfer and manufacture of digital prosthetic solutions such as single-tooth and bridge restorations on tiologic® TWINFIT implants.

CAD/CAM data sets for 3shape™, dental wings and exocad can be downloaded free of charge from www.dentaurem.com/digital. They include all the information needed for all types of restorations as a complete package, i.e. for printed models, hybrid abutments, one-piece abutments, bridges and bars.

5.1 Scan abutment titanium tiologic® TWINFIT, platform

The surface of the scan abutment titanium is optimized for digital capture, both intra-orally and on the model without scanning spray. The long cylinder at the connection point of the scan abutment titanium serves for easy integration into the (laboratory) implant and for exact axial detection of the position. The enclosed red retaining screw has a shortened thread, which only grips in the (laboratory) implant when the scan abutment has been inserted into the connection point in the correct position. The matching process of the scan data with the STL data takes place according to the specifications of the respective software manufacturer.

5.2. Laboratory implant tiologic® TWINFIT for printed and cast models
The laboratory implant can be used in both analog and digital workflow. The cavities of the laboratory analogs are stored in the CAD/CAM libraries of the various software manufacturers in order to achieve an optimum fit and the correct positioning of the laboratory analogs with the counter screw in 3D model printing. The parameters of the respective software must be checked and adjusted if necessary. In order for the counter screw supplied to fix the laboratory analog at the apical end in the printed model, the model base should be printed with a sufficient height. If a detachable mucosa is created, it should not be lower than 1.0 mm below the interface to ensure the optimal locking and retaining function of the laboratory analog in the cervical area of the printed model. The 3 retention tabs should still be visible in the printed model for monitoring purposes.

When the model is printed with the planned cavities, the laboratory implant for printed and cast models can be inserted in accordance with the intended series of abutments S, M or L. For this purpose, the laboratory implant is placed on the manual insertion key and fixed with the counter screw. The laboratory implant is then aligned in the cavity and pressed in. The closure screw supplied with the laboratory implant enables the laboratory implant to be fixed in the end position from apical.

5.3 CAD/CAM titanium block, tiologic® TWINFIT, PreForm

By using original Dentaurem CAD/CAM PreForm titanium blocks, it is possible to create highly precise abutments that are specific to each patient, have a precise fit and can be directly screw-retained. The original PreForm titanium blocks were designed for the Medentika® PreFace® abutment holders. These ensure a perfect connection to the tiologic® TWINFIT implant interface, both for in-house production and for production in a manufacturing center. CAD/CAM PreForm titanium blocks are available for all three series of abutments S, M and L with conical and platform connector geometry. Following selection of the respective indication in the data set in the corresponding CAD software, the design is made according to the specifications of the respective software manufacturer and according to the rules for dental prosthetics.

5.4 CAD/CAM titanium base tiologic® TWINFIT

The geometry of the titanium bases is designed to ensure a reliable, aesthetic bond with the ceramic mesostructure. tiologic® TWINFIT CAD/CAM titanium bases are available in all three series of abutments S, M and L with conical and platform connector geometry.

Following selection of the respective indication in the data set in the corresponding CAD software, the design is made according to the specifications of the respective software manufacturer and according to the rules for dental prosthetics.

The angulated screw aperture offers the optimum solution for the complete segment that is occlusally screw-retained. The screw aperture can be inclined at an angle of 20° to the implant axis, allowing discreet emergence in the palatal area, particularly in the anterior region.

CAUTION: Due to the higher gingival height in the conical design of the CAD/ CAM titanium base, a recess in the area to be bonded was constructed for the angulated screw aperture. This creates free space for a maximum angulation of 20°. In order to achieve the correct angle when planning, the scanning surface on the scan abutment must be opposite the desired exit of the angled screw aperture. The recess in the CAD/CAM titanium base (angulation alignment) is flush with a cam from the inner connection. In the case of the scan abutment, the upper scanning surface is located opposite a cam from the inner connection. During scanning, care must be taken to determine where the exit of the screw aperture should be located. For angulated screw apertures, the AnoTite screw for angulated screw apertures L 8.5 mm and the hexagonal ball end key 1.3 are used.

Before bonding, the head of the prosthetic screw, which fixes the CAD/CAM titanium base, is covered with wax. The bonding surfaces of the ceramic sleeve and of the CAD/CAM titanium base are blasted with aluminum oxide (50 µm/2 bar). It is also advisable to use the silanization method to condition this area on the CAD/CAM titanium base. Before bonding, ensure the surfaces are dry and free from grease. Adhesives must be used according to the manufacturer’s instructions (e.g.: "PANAVIA™ F2.0" by Kuraray Noritake or "Multilink® Hybrid Abutment" by Ivoclar Vivadent AG). Dentaurem GmbH & Co. KG gives no warranty for the stability of the meso/superstructure and the adhesive. After bonding, excess material is removed and the ceramic abutment is put into position (observing the rotational security). The all-ceramic crown is completed according to the manufacturer’s instructions.

Safety information.

- The wall thickness of the ceramic abutment must NOT be less than 0.5 mm.
- Preparation of a chamfer with angled inner edge and a minimum step of 0.5 mm.

5.5 CAD/CAM titanium base VARIO tiologic® TWINFIT

The titanium base VARIO has been specially developed for demanding cases in prosthetics.

In the case of an unfavorable implant position or in an aesthetically demanding area, there is both the possibility of laying the screw aperture towards oral and of shortening the shaft in various lengths according to the requirements. In the respective data set with the height of the adhesive cylinder, the adhesive surface to be milled in the ceramic abutment is automatically defined.

The titanium bases VARIO are available in all three series of abutments S, M and L with conical and platform connector geometry for tiologic® TWINFIT.

CAUTION: It is possible to open the adhesive cylinder on a pre-prepared surface for an angled screw aperture. When aligning the scan abutment and the CAD/CAM titanium base VARIO, it should then be noted that the alignment of the flat side on the sphere of the scan abutment and the exit direction of the angulated screw aperture (pre-prepared region) coincide.

Before bonding, the head of the prosthetic screw, which fixes the CAD/CAM titanium base, is covered with wax. The bonding surfaces of the ceramic sleeve and of the CAD/CAM titanium base are blasted with aluminum oxide (50 µm/2 bar). It is also advisable to use the silanization method to condition the surfaces to be bonded. Before bonding, ensure the surfaces are dry and free from grease. Adhesives must be used according to the manufacturer’s instructions (e.g.: "PANAVIA™ F2.0" by Kuraray Noritake or "Multilink® Hybrid Abutment" by Ivoclar Vivadent AG). Dentaurem GmbH & Co. KG gives no warranty for the stability of the meso/superstructure and the adhesive. After bonding, excess material is removed and the ceramic abutment is put into position (observing the rotational security). The all-ceramic crown is completed according to the manufacturer’s instructions.

Safety information.

- The wall thickness of the ceramic abutment must NOT be less than 0.5 mm.
- Preparation of a chamfer with angled inner edge and a minimum step of 0.5 mm.

5.6 CAD/CAM titanium adhesive bases

CAD/CAM titanium adhesive bases are suitable for adhesion in the case of multi-unit bridge restorations or full-arch restorations in the edentulous mandible or maxilla. They have a platform connector geometry for all series of abutments S, M and L which ensures the prosthesis sits optimally. There is no rotational security to enable maximum freedom when positioning the base on the implant. The cone for bonding has retention grooves. It is 3.0 mm high and can compensate for divergences of up to 30°. By using the bonding technique, the ceramic bridges, manufactured using CAD/CAM, can be fitted passively. Before bonding, the head of the prosthetic screw, which fixes the CAD/CAM titanium base, is covered with wax. The bonding surfaces of the ceramic sleeve and of the CAD/CAM titanium base are blasted with aluminum oxide (50 µm/2 bar). It is also advisable to use the silanization method to condition the surfaces to be bonded. Before bonding, ensure that the surfaces are dry and free from grease. Adhesives must be used according to the manufacturer’s instructions (e.g.: "PANAVIA™ F2.0" by Kuraray Noritake or "Multilink® Hybrid Abutment" by Ivoclar Vivadent AG). Dentaurem GmbH & Co. KG gives no warranty for the stability of the meso/superstructure and the adhesive. After bonding, excess material is removed and the ceramic abutment is put into position (observing the rotational security). The all-ceramic crown is completed according to the manufacturer’s instructions.

Safety information.

- NOT approved for single restorations.
- The wall thickness of the ceramic abutment must NOT be less than 0.5 mm.
- Preparation of a chamfer with angled inner edge and a minimum step of 0.5 mm.

Additional components

The prosthetic screw (not included in delivery) serves only to help the dental technician fix the abutment on the working cast. A new AnoTite screw should be used for the final fitting in the mouth.

Tightening torque

- Prosthetic screw, on model: manually
- AnoTite screw intra-orally: 30 Ncm
- AnoTite screw for angled screw aperture on the model: by hand
- AnoTite screw for angled screw aperture in the mouth: 25 Ncm

Attention: Do not confuse the prosthetic screw with the AnoTite screw. The prosthetic screw is colored silver and the AnoTite screw is anodized gold.

6. Composition

- CAD/CAM titanium base: titanium alloy (Ti6Al4V)
- CAD/CAM titanium adhesive base: titanium alloy (Ti6Al4V)
- CAD/CAM titanium block: titanium alloy (

- Hex keys:
 - Hexagon key - ISO shank 1.3: 1.4197
 - Hexagon key - ratchet 1.3: 1.4197
- Hexagonal ball end key - ISO shank 1.3: 1.4305**

Information about the composition is included in the material list for our implantology products (REF 989-801-05), or see www.dentaurum.com. The SSCP is available at https://ec.europa.eu/tools/eudamed.

7. Safety information

- The product should not be used if there is a known allergic reaction to one or more of the material components.
- Due to the small size, the article could be swallowed or aspirated. Aspiration could lead to difficulty in breathing or death due to asphyxiation. All articles used intra-orally should therefore be secured against swallowing and/or aspiration.
- All serious incidents arising from the use of the product should be reported to the manufacturer and the competent authority in the country in which the dental professional and/or the patient are resident.
- Different types of alloy in the oral cavity can lead to galvanic reactions.
- If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and must be disposed of if necessary.
- Only use tioLogic® TWINFIT prosthetic components in combination with tioLogic® TWINFIT implants.

8. Storage and shelf life

Store in a dry place.

9. Information for single use products

CAD/CAM titanium blocks, titanium bases (VARIO) and titanium adhesive bases are intended for single use only. Reconditioning of CAD/CAM titanium blocks, titanium bases (VARIO) and titanium adhesive bases (recycling) or reuse on patients are not permitted since it can then no longer be guaranteed that the article can be reprocessed safely or can function safely.

10. Information for multiple use products

The hex key 1.3, the hexagonal ball end key 1.3 and the scan abutment titanium are intended for multiple use. The product must be disinfected, cleaned and sterilized before being used again. Instruments and accessories must be replaced if they cannot be clearly identified or if the function is impaired due, for example, to poor readability of the markings or labels. Additional information can be found at www.dentaurum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

11. Information on disinfection/cleaning/sterilizing

tioLogic® TWINFIT prosthetic components are delivered non-sterile. They are for single use on one patient only. They must be cleaned and disinfected before being used on a patient. Dentaurum generally recommends that the product be sterilized in addition. Only sterilized prosthetic components may be used if bleeding occurs. Additional information can be found at www.dentaurum.com (Processing Instructions Prosthetic Components, Instruments and Accessories REF 989-801-09).

12. Product overview /scope of delivery, incl. REF

Please refer to the tioLogic® TWINFIT product catalog for the tioLogic® TWINFIT product range.

13. Quality information

Dentaurum ensures faultless quality of its products. The content of these Instructions for use is based upon our own experiences. The dental professional is solely responsible for the processing of the products. Liability for failures cannot be taken, as Dentaurum has no influence on the processing on site.

14. Explanation of symbols

i

Please refer to the label. Additional information can be found at www.dentaurum.com (Explanation of symbols REF 989-313-00).

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