## Estimado cliente

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## EN

### Dear customer,

Thank you for choosing a quality product from Dentaureum. Please read these Instructions for use to ensure you use this product in a safe and simple way and that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative.

#### 1. Manufacturer

Dentaureum GmbH & Co. KG  
Turnstr. 31 · 75228 Ispringen · Germany

#### 2. General device description

The hyrax® click screw (Hygienic Rapid Expansion Screw) is a special expansion screw for rapid palatal expansion. It has a screw locking feature (Fig. 1) to prevent unwanted reversing. The retention legs of this RPE screw are welded onto preformed bands (Fig. 2) or onto cap splints made of metal (Fig. 3). Alternatively they can also be polymerized into acrylic splints (Fig. 4) or anchored onto mini screws, e.g. tomas®-pin. The appliance is cemented or fixed to the tomas®-pin as usual. Follow the instructions from the cement manufacturer.

#### 3. Intended purpose

Used for rapid palatal expansion with a fixed appliance.

#### 4. Intended users

All products are intended for dental professionals only.

#### 5. Indications

In case of transverse under-development of the upper jaw.

#### 6. Contraindications

Stainless steel expansion screws must not be used in patients with a nickel and/or chromium allergy. Warnings are found on the respective product labels. Also see 9. Composition.

#### 7. Patient target group

All patients for whom a fixed appliance with an expansion screw can be placed as part of orthodontic treatment.

#### 8.1. Fabrication of hyrax® click appliance with preformed bands

Preformed bands are fitted on the upper jaw left and right first bicuspid and left and right first molars. Using an impression, these bands are either transferred to or fitted on the hard plaster model. The retention legs of the hyrax® click screw are then bent so that they can be welded onto the four adapted bands (Figs. 2 + 5).

The hyrax® click screw should be "suspended" at least 1 mm – 2 mm above the palatal vault/gingiva (cf. Fig. 6) and its dorsal side should be in line with the mesio-palatal cusps of the first molars. In order to avoid leverage, the screws must be situated parallel to the occlusal plane. The hyrax® click screw is positioned so that it may be turned in a dorsal direction (Figs. 2 – 4) in order to make the task of inserting the safety key and activating the screw easier.

#### 8.2. Tips for working with hyrax® click screws

The bands can be connected by means of a remalay® wire ø 1.0 mm (REF 528-100-00) on the palatal and the buccal side as reinforcement. These reinforcement wires are laser-welded above the retention legs. The palatal reinforcement wires should extend slightly over the bands and be supported by the adjacent teeth in order to distribute the expansion forces more effectively.

In order to bend the retention legs of the hyrax® click screws gently and easily, we recommend using the universal bending key (REF 611-122-00, Figs. 7 + 9) without the aid of additional pliers (lever principle).

If this key is not available, the retention legs must be stabilized by holding the expansion screw body with a pair of sturdy pliers, e.g. wire bending pliers Maxi Premium-Line (REF 013-522-00), in order to prevent the retention legs from breaking (Fig. 8).

#### 9. Composition

The hyrax® click screws are made of the following materials: DIN 1.4301, DIN 1.4305 and DIN 1.4310. The composition is included in the materials list; please refer to the catalog or [www.dentaureum.com](http://www.dentaureum.com). The SSCP is available at <https://ec.europa.eu/tools/eudamed> and [www.dentaureum.com](http://www.dentaureum.com).

#### 10. Safety information

Always use only the blue safety key supplied to prevent swallowing (Fig. 10). The product should not be used if there is a known allergic reaction to one or more of the material components.

If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and if necessary must be disposed of.

#### 11. Information for single use products

The product is intended for single use only. Reconditioning (recycling) of the product and its reuse is not permitted, as perfect functionality of the product cannot be guaranteed.

#### 12. Further information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident are informed accordingly.

#### 13. Notes on disposal

When disposing of the product, please adhere to the national regulations that are valid for your region and - if available - the details outlined in the safety data sheet

#### 14. Product range

For the complete range of hyrax® click screws, please see our orthodontics catalog or visit [www.dentaureum.com](http://www.dentaureum.com).

#### 15. Quality

Dentaureum ensures faultless quality of its products. These recommendations are based upon Dentaureum's own experiences. The dental professional is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as Dentaureum has no influence on the processing on site.

#### 16. Explanation of symbols used on the label

⚠ Please refer to the label. These Instructions for use and additional information can be found at [www.dentaureum.com](http://www.dentaureum.com) (Explanation of symbols REF 989-313-00).

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