Estimado cliente

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EN

Dear customer,

Thank you for choosing a quality product from Dentaurum. Please read these instructions for use to ensure you use this product in a safe and simple way and that you and your patients gain full benefit. In case of questions or ideas, please contact your local representative.

1. Manufacturer

Dentaurum GmbH & Co. KG
Turnstr. 31 | 75228 Ispringen I Germany

2. General device description

hyrax® and hyrax® neo screws without retention legs (hygienic rapid expansion screw, Figs. 1–2) are special screws used for fast palatal expansion. The body of the RPE screw is welded onto individual retention legs, connectors or cap splints which have been either cast or produced by the SLM method in metal. The appliances are cemented or fixed to the tomas®-pin as usual. Follow the instructions from the manufacturer.

3. Intended purpose

Used for fast palatal expansion with a fixed appliance.

4. Intended users

All products are intended for dental professionals only.

5. Indications

In case of transverse under-development of the upper jaw.

6. Contraindications

Stainless steel expansion screws must not be used in patients with a nickel and/or chromium allergy. Warnings are found on the respective product labels. Also see 9. Composition.

7. Patient target group

All patients for whom a fixed appliance with an expansion screw can be placed as part of orthodontic treatment.

8.1. Fabrication of a hyrax®/hyrax® neo appliance with individualized retention legs and cap splints

The following steps should be followed to fabricate individualized retention legs and cap splints by casting or by the SLM method (Figs. 3 – 5). Individualized cap splints are fitted on the first premolars in the maxilla and the first molars, each left and right (Figs. 10 – 11). The retention legs are connected to the hyrax®/hyrax® neo screw in such a way that two weld seams form along the lateral or the upper groove (Figs. 8 – 9). The hyrax®/hyrax® neo screw should be “suspended” at least 1 mm – 2 mm above the palatal vault/gingiva (Figs. 6 – 7) and its dorsal side should be in line with the meso-palatal cusps of the first molars. In order to avoid leverage, the screws must be situated parallel to the occlusal plane. The hyrax®/hyrax® neo screw is positioned so that it may be turned in a dorsal direction in order to make the task of inserting the safety key (Fig. 13) and activating the screw easier.

8.2. Tips for working with hyrax®/hyrax® neo screws without retention legs

For the connection to be durable, the individualized retention legs must be welded onto the groove of the hyrax®/hyrax® neo screw body (Figs. 8 – 9), for example using the desktop Compact welding unit (REF 090-578-50). As an average value for NiCr, the initial laser setting on the unit is as follows: 240 volt, 4 ms, 0,8 mm. The use of a shielding gas is indispensable during the welding process. When using or adding new material to reinforce the welding area, we recommend the NiCr welding wire, ø 0,50 mm (REF 528-220-00, Fig. 12).

9. Composition

The hyrax® and hyrax® neo screws are made of the following materials: DIN 1.4301, DIN 1.4305 and DIN 1.4310. The material composition is included in the orthodontics materials list; please refer to the orthodontics catalog or www.dentaurum.com. The SSCP is available at https://ec.europa.eu/tools/eudamed and www.dentaurum.com.

10. Safety information

Always use only the blue safety key (Fig. 13) supplied to prevent swallowing!

The product should not be used if there is a known allergic reaction to one or more of the material components.

If the packaging shows serious damage, the product must be checked for integrity and cleanliness before use, and if necessary must be disposed of.

11. Information for single use products

The product is intended for single use only. Reconditioning (recycling) of the product and its reuse is not permitted, as perfect functionality of the product cannot be guaranteed.

12. Further information

Should the dental professional and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the dental professional and/or the patient is resident are informed accordingly.

13. Notes on disposal

When disposing of the product, please adhere to the national regulations that are valid for your region and - if available - the details outlined in the safety data sheets.

14. Product range

For the complete range of hyrax®/hyrax® neo screws, please see our orthodontics catalog or visit www.dentaurum.com.

15. Quality information