
Gebrauchsanweisung | Instructions for use |
Mode d'emploi | Modo de empleo

Paratherm

REF 094-270-00

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Dear customer

Thank you for choosing a quality product from Dentaureum.

It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

Instructions for use cannot describe every eventuality and possible application. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaureum.com, even though you frequently use the same product.

1. Manufacturer

Dentaureum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

2. Safety instructions

 **Caution** **Read through the instructions carefully before installing and operating the unit. Only then should the unit be switched on.**

2.1 Disposal

Disposal regulations of the country of operation apply when a Dentaureum product is taken out of operation once and for all. Dentaureum or the dental trade is available to answer questions regarding the correct disposal of any specific product.

Note:

This symbol indicates that the products should not be disposed of together with household garbage. The legislator disallows commercial customers to return electronic waste via municipal collection points. Further information is available from Dentaureum or the dental trade.

2.2 Delivery

| | Description | | Quantity |
|---|---|--|----------|
| | Paratherm unit with surveying / blocking-out arm | | 1 |
| Ⓔ | Thermal blocking-out instruments 0°, 2° | | 1 each |
| Ⓓ | Undercut gauges 0.25 / 0.35 / 0.50 / 0.75 mm | | 4 |
| Ⓐ | Analyzing rod | | 1 |
| Ⓑ | Graphite holder | | 1 |
| Ⓒ | Replacement graphite markers | | 10 |
| | Plastic cover | | 1 |
| | Fuse T 0.25 A | | 1 |
| | Instructions for use with EEC-Declaration of conformity | | 1 |

3. Area of use

Cast partial denture, combination techniques, conical crowns, positioning of snap attachments and attachments.

4. Installation


Set the necessary voltage on the underside of the unit (110 / 220 / 240 V). Then connect the unit to the mains supply. Stand the Paratherm on a level working surface.

5. Operating the Paratherm

Mount the master model on the model table already available or delivered with the unit. Loosen the knurled screw ② and raise the sliding rods to the necessary level. Fasten with retainer ring ①. Clamp the instrument, e.g. analyzing rod, in the surveying handpiece ④ and push the selected blocking-out instrument into the thermal blocking-out handpiece ⑤.

Place the model table ⑥ in the most convenient position on the work surface ⑦. Loosen the locking screws ③ so that the arms may be moved freely. Switch on the thermal blocking-out handpiece ⑧ and set the desired temperature (maximum 160 °C) on the controller ⑨.

! Caution: Only change the blocking-out instruments ⑩ whilst they are still warm.

5.1 Electric wax knife

The knife is connected by means of the socket on the rear side of the unit. Only operate the wax knife with an inserted knife. Turn the power switch ⑧ on.

Set the desired temperature on the potentiometer ⑪. Maximum temperature 280 °C.

! Caution! Only change the knife whilst it is still warm.

5.2 Locking the surveying arm

Rotate the surveying handpiece ④ forwards of the centre of the work surface and lock into position with screws ②, ③.

Clamp the wire aid in the surveying arm and fix the level using ⑫. With the model already clamped into position, guide the model table onto the wire aid. Correct the level.

When ordering spare parts please quote the unit number (type plate on rear side of the unit).

6. Spare parts

| | Description | | Quantity | REF | |
|---|---------------------------------|---------|----------|------------|--|
| Ⓔ | Thermal blocking-out instrument | 0° | 1 | 094-275-00 | |
| | | 2° | | 094-277-00 | |
| Ⓓ | Undercut gauge | 0.25 mm | 1 | 094-283-00 | |
| | | 0.35 mm | | 094-284-00 | |
| | | 0.50 mm | | 094-285-00 | |
| | | 0.75 mm | | 094-286-00 | |
| Ⓐ | Analyzing rod | | 1 | 094-287-00 | |
| Ⓑ | Graphite holder | | 1 | 094-288-00 | |
| Ⓒ | Replacement graphite markers | | 10 | 094-289-00 | |

7. Recommended accessories

(Not included in unit REF 094-270-00)

| | Description | | Quantity | REF | | |
|---|--|----|----------|------------|--|--|
| Ⓔ | Thermal blocking-out instrument | 1° | 1 | 094-276-00 | | |
| | | 3° | | 094-278-00 | | |
| | | 4° | | 094-279-00 | | |
| | | 6° | | 094-280-00 | | |
| | | 8° | | 094-281-00 | | |
| | Model table | | | 086-340-00 | | |
| | <p>Electric wax knife Handpiece, cable and 1 wax knife, holder for handpiece.</p> <p>In order to regulate the temperature, a transformer and potentiometer have been built into the Paratherm (REF 094-270-00)</p> | | | 094-272-00 | | |

| | | | | | |
|--|--------------------------------|--|--|------------|---|
| | Modelling knife, small | | | 089-602-00 |  |
| | Thomas modelling point, fine | | | 089-603-00 |  |
| | Thomas modelling point, medium | | | 089-604-00 |  |
| | Orthodontic wax knife | | | 089-605-00 |  |
| | Modelling spatula | | | 089-606-00 |  |

8. Quality information

Dentaurum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user is responsible for the processing of the products. Responsibility for failures cannot be taken, as we, Dentaurum, have no influence on the processing on site.

9.

EC-Declaration of Conformity

DENTAURUM GmbH & Co. KG
Turnstr. 31
75228 Ispringen

hereby declares that the design and construction of the laboratory equipment described below, including the version marketed by us, comply with the basic regulations governing safety and health as stated in the EC Guidelines. This declaration will become invalid if the laboratory equipment is modified or altered in any way without our prior consent.

| | | |
|----------------------------|----------------------|---|
| Description of unit: | Paratherm | |
| REF: | 094-270-00 | |
| Start with unit No.: | 32-080 | |
| EC guidelines: | 89/336/EEC | Electromagnetic compatibility |
| | 92/31/EEC | Modified Guideline |
| | 73/23/EEC | Electrical equipment used within certain voltage limits |
| | 93/68/EEC guidelines | Modifications of the above |
| Applied unified standards: | EN 55011 | |
| | EN 50082 | |
| | EN 61010 | |

Date and manufacturers signature: 01.01.2010
Signatory:

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