remanium® 

CE 0483

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Dear customer,

Thank you for choosing a quality product from Dentaureum.

It is essential to read these Processing Instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit.

It is not possible to describe every eventuality and possible application in these instructions. In case of questions or ideas, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaureum.com, even though you may frequently use the same product.

Manufacturer

Dentaureum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

Quality

Dentaureum ensures faultless quality of its products. These recommendations are based upon Dentaureum's own experiences. The user is solely responsible for the processing of the products. Responsibility for failures cannot be taken, as we, Dentaureum, have no influence on the processing on site.

Contents

1. General instructions for processing	14
2. Wax up	14
3. Spruing	14-15
4. Investing	14
5. Casting ring systems	16
6. Wax burn-out and preheating	16
7. Melting and casting	16
7.1. Crucible	16
7.2. Metal quantity	16
7.3. Reuse of casting buttons	16
7.4. Cast-on	17
7.5. Casting equipment	17
8. High-frequency melting with vacuum pressure casting units	17
8.1. Premelting	17
8.2. Melting	17
9. High-frequency melting with centrifugal casting units	17
10. Melting with torch	18
10.1. Observing the melt during torch melting	18
11. Overheating	18
12. Cooling of the casting ring	18
13. Devesting and sandblasting	19
14. Finishing and polishing	19
14.1. Finishing and polishing steps	19
14.2. Electrolytic polishing	20
15. Laser welding	20
16. Soldering	20
17. Cleaning	21
18. Contraindications and adverse reactions	21

1. General instructions for processing

These Processing Instructions deal with the main steps when processing non-precious remanium® cast partial denture alloys.

Individual additional information regarding our non-precious remanium® cast partial denture alloys which may differ from these instructions can be found in the information supplied in each individual alloy package or in the internet under www.dentaurum.com.

In case of further questions, please contact your local representative.

Please refer to the information at the end of these instructions for information on contra-indications and adverse reactions.

2. Wax-up

The outstanding properties of non-precious remanium® cast partial denture alloys allow for very thin frameworks. A thickness of 0.4 mm is sufficient for full metal dentures. If the transversal connectors on partial dentures are slender in design, then the main body structure must be stronger. A cross-section diameter of 2 x 4 mm is recommended for sublingual bars.

Waxes and acrylics which burn out without leaving a residue are essential for high quality casting results.

3. Spruing

Round profiles \varnothing 3.5 – 4.0 mm. Length 20 mm minimum. For voluminous wax ups \varnothing 4.0 mm. For upper jaw 2 – 4 sprues, for lower jaw 2 – 3 sprues.

The sprues should be attached to the pattern in such a way as to determine the flow direction of the molten metal. It is possible to cast from above as well as through the model itself.

The overview on page 15 shows recommendations of sprue diameter and position depending on different types of design and construction.

4. Investing

Use phosphate-bonded investment materials with optimally adapted expansion, such as rema® Exakt, rema® Exakt F or rema® dynamic S. rema® dynamic S is an investment material recommended for quick casting methods with which it is possible to immediately place the ring into the furnace at final temperature without risking a bad fit. For 3D printed objects, we recommend rema® CAD/Vest.

Spruing

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

3.5 mm

4 mm

4 mm

2.5 mm / 4 mm

4 mm

3.5 mm

3.5 mm

4 mm


2.5 mm / 4 mm

2.5 mm / 4 mm


2.5 mm / 4 mm

3.5 mm

25

3.5 mm

26

Length min. = 20 mm

Length max. = 40 mm

5. Casting ring systems

It is recommended to use reusable casting rings when investing, for example rema® Form quick investment system or Neo-Star investment system.

The model base must always be surrounded by investment material. This prevents the model from being separated from the ring during the casting process.

6. Wax burn-out and preheating

Observe the Instructions for use for the investment material used. The ring temperature of the investment material should be kept as low as possible. The temperatures range between 850 °C / 1562 °F and 1050 °C / 1922 °F according to the type of investment, casting machine and pattern design used.

A lower preheating temperature provides a more homogenous metal structure and reduces oxidation. When deciding on a preheating temperature it is also important to take into account the pressing pressure of the casting machine, the thickness of the wax pattern and the investment material's gas permeability.

7. Melting and casting

7.1. Crucible

Use a ceramic-based crucible (magnesium oxide, silicon oxide or aluminum oxide) that is suitable for non-precious alloys.

Do not use graphite crucibles.

Preheat empty crucible in the furnace. Use one crucible per alloy and clean after each casting.

Do not use special flux powder!

For larger quantities of metal, a preliminary melting phase (casting delay phase) is recommended.

7.2. Metal quantity

In general, one cast partial denture requires 4 – 5 casting ingots at 6 g each.

7.3. Reuse of casting buttons

The alloy may be changed in its composition and properties if reused. Its reuse is therefore not permitted. This is the only way to ensure clear batch traceability.

7.4. Cast-on

Metals used for cast-on and high platinum content attachments may not be used in combination with non-precious remanium® cast partial denture alloys.

7.5. Casting equipment

Melting and casting techniques suitable for casting non-precious remanium® cast partial denture alloys include high-frequency melting with vacuum pressure or centrifugal casting units, arc melting with vacuum pressure casting, or torch melting with propane/oxygen or acetylene/oxygen or other gas mixtures capable of reaching the casting temperature.

8. High-frequency melting with vacuum pressure casting units

Set vacuum to 250 – 450 mbar (recommended pressure: 450 mbar).

8.1. Premelting

Premelt the ingots until they are red hot (no visible rounding at the corners of the ingots).

8.2. Melting

Depending on the type of alloy, melt the metal until the skin which forms over the molten metal breaks open, **then cast immediately** or wait for 1 – 2 seconds after the oxide skin layer has broken open and then cast (see the relevant alloy Instructions for use).

Observe the melting process and initiate the casting process manually. Oxides formed during melting (see melting illustrations in the operating instructions for the pressure casting machine) remain as slag in the crucible after casting and can be easily removed. If the correct melting and casting technique is used, there is only a minimum residue in the crucible and the molten metal is not overheated.

Automatic casting with a pyrometer control is unsuitable, as the pyrometer cannot always accurately determine the right time for casting.

9. High-frequency melting with centrifugal casting units


Once all ingots have melted completely, wait until the skin which forms over the molten metal breaks open, **then cast immediately** or wait for 1 – 2 seconds after the oxide skin layer has broken open and then cast (see the relevant alloy Instructions for use).

We recommend premelting the metal if a large amount is to be cast.

10. Melting with torch

When melting, the remanium® cast partial denture alloys show a clean melt without sparks.

Caution! Use torch with shower head nozzle!

a) Propane / oxygen melting

Select maximum torch setting as stated by manufacturer.
Oxygen 2 – 3 bar / propane: 1 bar.

b) Acetylene/oxygen melting

Follow the instructions of the torch manufacturer.
Acetylene: 0.7 bar / oxygen: 3 bar.

The blue core of the flame at the torch head should be 4 – 5 mm in length. Preheat the empty crucible. Melt the metal in the crucible with a light circular movement of the flame. The distance at which the torch nozzle should be held depends upon the type of nozzle and the gas mixture used.

10.1. Observing the melt during torch melting

When melting with a torch, an oxide layer forms when the last piece of metal has collapsed. With a light circular movement of the flame, continue the melting process until the metal under the oxide skin is seen to move under the pressure of the flame. Begin the casting process **before** the oxide layer breaks open.

11. Overheating

As with all melting and casting techniques, **overheating the metal must be avoided**. Overheating results in blowholes, micro-porosities and grain coarsening. These casting imperfections are frequently reasons for cracks.

12. Cooling of the casting ring

Remove the casting ring carefully and allow to cool. In order to achieve the best mechanical qualities and a stress-free material allow the metal to cool down slowly to room temperature.

After at least 30 minutes of bench cooling, the casting ring may, if necessary, be carefully quenched in water.

Caution! Danger of fire. Always place hot casting rings onto heat-resistant surfaces.

Ensure that a safe distance from flammable objects is observed.

13. Devesting and sandblasting

Avoid dust inhalation!

Please observe the recommendations made on safety precautions in the Instructions for use enclosed with the product. Moisten the casting ring before devesting in order to reduce the amount of dust produced. Do not strike the casting button with a hammer. Sandblast using pure aluminum oxide, grain size approx. 125 µm to 250 µm.

14. Finishing and polishing

14.1. Finishing and polishing steps

Steps	Tools	REF	Comments
Removal of casting sprues	Supercut STM separating disc	130-102-00	High cutting efficiency and long life span for fast cutters.
	Supercut STM separating disc	130-111-00	High cutting efficiency and long life span for handpiece.
	Supercut STM separating disc	130-113-00	High cutting efficiency and long life span for handpiece.
Grinding of sprues	Grinding wheel	131-322-50	Highly suitable for grinding the stump of the sprue and for coarse grinding.
Finishing	Hard metal bur, with cross cut	123-581-00	General use in inaccessible areas.
		123-583-00	
		123-584-00	
123-585-00			
	Diasint burs, instruments	137-826-00	Sintered diamond burs. Metallic bond prevents fragmentation.
		137-827-00	
		137-831-00	
	Grinding points	132-802-00	Ceramic-bonded corundum grinding points.
		132-803-00	
Electrolytic polishing	Electropolishing solution	128-301-00	Cover attachments.
Buffing	Silichrom polisher	138-640-00	Abrasive performance, rough.
	Green rubber polishers	138-645-00	Abrasive effect, fine.
		138-101-00	
Polishing	Polishing brushes	141-712-00	Universal polish with hand tool.
	Tiger Brillant polishing paste	190-350-00	Fast acting high shine polishing paste.
	Tiger Starshine Universal finishing paste	190-301-00	High-polishing finishing paste.

General information:

Please observe the recommendations made on safety precautions to be taken when processing the alloy mechanically. These can be found in the Instructions for use enclosed with the product.

Dentaurum's remanium® cast partial denture alloys can be finished easily and safely with special finishing and polishing instruments. Suitable instruments and materials are available in the Dentaurum product range.

Successful fitting and finishing of the remanium® cast partial denture alloys depend greatly upon the accuracy of the wax pattern and casting using Dentaurum investment materials which have been optimally adapted for these alloys.

Use the specific instruments for only one alloy.

14.2. Electrolytic polishing

The duration of the polishing times depends upon the type of polishing unit and electrolyte used. General duration: Polishing time without movement in the solution approx. 2 x 5 minutes (after half time change the position of the object in the solution) at 4 amps, and approx. 10 minutes at 4 amps with movement.

Polishing baths should be warmed in order to increase their polishing ability (30 °C / 86 °F to 40 °C / 104 °F). Hot baths of > 50 °C / 122 °F become damaged and diminish the polishing effect.

15. Laser welding

All remanium® alloys are suitable for laser welding due to their extremely low or non-existent carbon content. remanium® GM 900, free of carbon, is especially suitable.

The laser welding technique produces a solder-free, mechanically strong and corrosion-resistant joint.

It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters of each individual laser machine. Suitable welding wires are available as filler materials for all remanium® CoCr alloys.

CoCr welding wire:

ø 0.25 mm	Roll 2 m long	REF 528-215-10
ø 0.35 mm	Roll 2 m long	REF 528-210-10
ø 0.5 mm	Roll 2 m long	REF 528-200-10

16. Soldering

Try to avoid soldered joints if possible.

If however, it is necessary to solder, please ensure a suitable solder with the correct composition and melting temperature which is appropriate for the alloy is used.

Coat the solder objects with an adequate amount of flux!

Caution! Avoid letting the flux dry out!

Heat the framework until it is red hot. Position the parts that have been dipped in the flux and heat to flowing.

The following are suitable solders and flux for remanium® CoCr alloys:

Solder	CoCrMo-Sold 1	approx. 3 g	REF 102-306-00
Flux	rema® Flux 1	25 g	REF 102-304-00

17. Cleaning

remanium® frameworks can be cleaned with commercially available cleaning agents for dentures. Refer to the cleaning recommendations on the labels.

The use of household cleaning agents is not recommended as these could cause the frameworks to discolor.

Recommended ultrasonic solution after polishing:

Ultrasonic cleaning solution	1000ml	REF 094-610-00
------------------------------	--------	-----------------------

18. Contraindications and adverse reactions

Signs of intolerance to non-precious remanium® cast partial denture alloys are extremely rare if the manufacturer's Instructions for use are adhered to.

If the patient has a proven allergy against any component within the alloy, this alloy must not be used for safety reasons.

There have been individual reportings of local irritations which were electrochemically induced.

If various alloys have been used, it is possible that galvanic effects may occur.

Please observe the information contained in the enclosed Instructions for use.

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