



discovery[®]
delight



Surprisingly different.
The lingual system.

D
DENTAURUM


Incredibly effective.



discovery® delight dramatically reduces the number of components required for treatment while ensuring a free selection of all treatment options, from 2D to 3D treatment.

In collaboration with Prof. Dr. Rubens Demicheri, Montevideo, a lingual system that addresses patient needs and offers all the advantages of the well-proven discovery® bracket family emerged:

- **highest quality**
- **virtually invisible**
- **short treatment time**
- **high wearing comfort**
- **perfect design**




Simply practical.



The main components of **discovery® delight** are four bracket types: they can be used both for the upper and lower jaw and ensure easy and time-saving handling!

The brilliant combination of a vertical slot in the anterior region with a horizontal slot in the posterior region (open / closed) enables optimal force transfer.

Carefully selected arch wires and instruments round off the product portfolio of **discovery® delight**.



Anterior bracket

Posterior bracket
open



1-piece bracket

High biocompatibility

Laser-structured base

Optimal retention

Optimal wing design

Easy and safe ligating

Delicate and low profile

Largest possible space for the tongue

Posterior bracket
closed

Molar tube

Remarkably versatile.



Findings

Minor position deviation
of the posterior teeth
(social six)

Treatment

2D

**Average
treatment time**

approx. 6 months



Moderate overcrowding in the anterior region without rotation of the premolars

Severe overcrowding in the anterior region without rotation of the premolars

Severe overcrowding in the anterior region with rotation of the premolars

Position deviations in the entire dental arch

2D

2D

3D

3D


approx. 12 months

approx. 6 – 24 months

approx. 6 – 24 months

approx. 6 – 24 months

Clearly arranged.



*The whole is greater
than the sum of its parts.*

Aristoteles



Lingual brackets /
lingual tubes

- Anterior bracket
3-3/Upper – Lower
- Posterior bracket , open
4-5/Upper – Lower
- Posterior bracket , closed
4-5/Upper – Lower
- Molar tube
6er/Upper – Lower



Lingual arches round

- rematitan® LITE, ø 10
- rematitan® LITE, ø 12
- rematitan® LITE, ø 14
- rematitan® LITE, ø 16
- rematitan® LITE, ø 18
- remanium®, ø 16
- remanium®, ø 18
- dentaflex®, ø 15

Lingual arches rectangular

- rematitan® LITE, 16 x 16
- rematitan® LITE, 19 x 19
- remanium®, 18 x 18
- remanium®, 19 x 19



Instruments

- Positioning tool
- Mirror with auxiliary lines
- Ligature cutter
- Distal end cutter
- Crimping tweezers

Rich selection of brackets.



In addition to the new **discovery® delight** system, you can make your selection among other brackets of the discovery® family. They meet all treatment needs and offer numerous combination options within this range.

Brackets for any treatment!



discovery®
pearl

Fascinatingly natural.



discovery®
smart

Small and brilliant.



discovery[®]

The sum of experience.



discovery[®] *sl* *discovery*[®]
sl 2.0

Proven technology for your success.

Dentaurum Group

Germany | Benelux | España | France | Italia | Switzerland | Australia | Canada | USA
and in more than 130 countries worldwide



DENTAURUM
QUALITY
WORLDWIDE
UNIQUE



All clinical figures provided by courtesy of Prof. Rubens Demicheri, Dr. med. sci

Date of information: 05/19
Subject to modifications

D
DENTAURUM

Turnstr. 31 | 75228 Ispringen | Germany | Phone +49 72 31/803-0 | Fax +49 72 31/803-295
www.dentaurum.com | info@dentaurum.com