



Assortments. Lingual brackets.

Case	REF	Contents	Brackets*	REF
1. Social six pack	891-406-00	6 x	Anterior bracket	790-450-00
2. Single jaw	891-407-00	6 x	Anterior bracket	790-450-00
		4 x	Posterior bracket	790-452-00
3. Single jaw + tubes	891-408-00	6 x	Anterior bracket	790-450-00
		4 x	Posterior bracket	790-452-00
		2 x	Molar tube	724-102-51
4. Complete case	891-409-00	12 x	Anterior bracket	790-450-00
		8 x	Posterior bracket	790-452-00
		4 x	Molar tube	724-102-51

* Slot: \square 0,51 x 0,51 mm / 20 x 20



The arch sequence.

Lingual arches.

Case	Seq.	Diameter	Material	REF
1. Social six	1.	∅ 0,38 mm / 15	dentaflex® lingual arch, round, 6-strand twisted	766-771-00
	2.	∅ 0,25 mm / 10	rematitan® "LITE" lingual arch, round	766-760-00
	3.	∅ 0,30 mm / 12	rematitan® "LITE" lingual arch, round	766-761-00
	4.	∅ 0,40 mm / 16	remanium® lingual arch, round	766-767-00
2. Single jaw	1.	∅ 0,38 mm / 15	dentaflex® lingual arch, round, 6-strand twisted	766-771-00
	2.	∅ 0,25 mm / 10	rematitan® "LITE" lingual arch, round	766-760-00
	3.	∅ 0,30 mm / 12	rematitan® "LITE" lingual arch, round	766-761-00
	4.	∅ 0,35 mm / 14	rematitan® "LITE" lingual arch, round	766-762-00
	5.	∅ 0,40 mm / 16	rematitan® "LITE" lingual arch, round	766-763-00
	6.	∅ 0,40 mm / 16	remanium® lingual arch, round	766-767-00
3. Single jaw + tubes	1.	∅ 0,38 mm / 15	dentaflex® lingual arch, round, 6-strand twisted	766-771-00
	2.	∅ 0,25 mm / 10	rematitan® "LITE" lingual arch, round	766-760-00
	3.	∅ 0,30 mm / 12	rematitan® "LITE" lingual arch, round	766-761-00
	4.	∅ 0,35 mm / 14	rematitan® "LITE" lingual arch, round	766-762-00
	5.	∅ 0,40 mm / 16	rematitan® "LITE" lingual arch, round	766-763-00
	6.	∅ 0,41 x 0,41 mm / 16 x 16	rematitan® "LITE" lingual arch, rectangular	766-765-00
	7.	∅ 0,46 x 0,46 mm / 18 x 18	remanium® lingual arch, rectangular	766-769-00
4. Complete case	1.	∅ 0,38 mm / 15	dentaflex® lingual arch, round, 6-strand twisted	766-771-00
	2.	∅ 0,25 mm / 10	rematitan® "LITE" lingual arch, round	766-760-00
	3.	∅ 0,30 mm / 12	rematitan® "LITE" lingual arch, round	766-761-00
	4.	∅ 0,35 mm / 14	rematitan® "LITE" lingual arch, round	766-762-00
	5.	∅ 0,40 mm / 16	rematitan® "LITE" lingual arch, round	766-763-00
	6.	∅ 0,41 x 0,41 mm / 16 x 16	rematitan® "LITE" lingual arch, rectangular	766-765-00
	7.	∅ 0,48 x 0,48 mm / 19 x 19	rematitan® "LITE" lingual arch, rectangular	766-766-00
	8.	∅ 0,48 x 0,48 mm / 19 x 19	remanium® lingual arch, rectangular	766-770-00