

tiologic
TWINFIT

conical platform

IT'S MY CHOICE.

ABUTMENT SWITCH

2 CONNECTOR GEOMETRIES – 1 IMPLANT.

conical platform





Depth-stop system.

Flexible, safe preparation.

With the depth-stop system for situational single use, surgical preparation of the site is safe. Depth-stop sleeves, if required, can be pushed onto the drill quickly and easily. Color coding states the implant diameter. The sleeves are sterile in the packaging and are included with every implant.

- Safe implant insertion is possible with or without depth-stop
- Depth-stop sleeves easy to incorporate
- Preparation gentle on the bone
- Color coding for diameters



Packaging.

Easy, efficient and safe to handle.

- PickUp for direct, safe grasping of implant
- Color-coded depth-stop sleeves for single use included
- Closure screw included
- Double sterile packaging



The new, patented implant system tiologic® TWINFIT with maximum flexibility from insertion to prosthetic restoration. Flexible and individual; choose the type of abutment best for your patient – conical or platform.

- Conical or platform – Abutment Switch
- 2 connector geometries on one implant
- Depth-stop system for situational single use
- Comprehensive assortment of prosthetics to suit all requirements
- CAD/CAM, optimized for a flexible, digitized workflow




Angulated solutions for all requirements.

4Base system.

- Screw-retained superstructures
- Uniform interface facilitates prosthetic restoration
- Angulations up to 50° possible
- Avoids augmentative measures, shortens time in chair

Angulated screw aperture.

- Counteracts implant divergencies
- For aesthetic, functional prosthetic results
- Discreet access in anterior region
- One screwdriver for all angulated screw apertures



Digital workflow.

Complete solution for your CAD/CAM process.

The patented implant system tiologic® TWINFIT is designed to suit a digital workflow. As your digital partner, the Dentaurem Group places great value on process sequences that are efficient and easy to understand, using materials from us – from scanning through to manufacture. Sequences can, of course, be adapted to suit your individual needs.



Dentaurum Group

Germany | Benelux | España | France | Italia | Switzerland | Australia | Canada | USA
and in more than 130 countries worldwide.



DENTAURUM
QUALITY
WORLDWIDE
UNIQUE

➔ For more information on our products and services, please visit www.dentaurum.de

Date of information: 11/18

Subject to modifications



www.dentaurum.com



Like us on Facebook!



YouTube

Visit us on YouTube!



Follow us on Pinterest!

